RUSH ROGUE® 2 by PROTEOR

ROG2 / EVQR2 / H2R2

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EN Instruction for use

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NO Brukerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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	RUSH ROGUE® 2 by PROTEOR <i>Instructions for use</i> Read before use	IFU-01-101 Rev. D 2022-11
---	--	---------------------------------

Provide patient with safety instructions and information necessary to use the device (§11, 12, 13, 14, 15, 17, 18, 19)

1. INCLUDED ITEMS

RUSH ROGUE 2 FOOT		
Part description	Part number	Included / Sold separately
RUSH ROGUE 2 Foot with regular toe	ROG2-xx-y-RU	Included
RUSH ROGUE 2 Foot with sandal toe, right	ROG2-xx-y-SR	
RUSH ROGUE 2 Foot with sandal toe, left	ROG2-xx-y-SL	
RUSH Standard Foot Shell, With CAP	FS	Sold separately (see catalogue)
RUSH Sandal Toe Foot Shell, No CAP	FTC-4R-1	
Black Spectra® sock	SS	Included
RUSH ROGUE 2 EVAQ8		
Part description	Part number	Included / Sold separately
RUSH ROGUE 2 EVAQ8 Foot with regular toe	EVQR2-xx-y-RU	Included
RUSH ROGUE 2 EVAQ8 Foot with sandal toe, right	EVQR2-xx-y-SR	
RUSH ROGUE 2 EVAQ8 Foot with sandal toe, left	EVQR2-xx-y-SL	
RUSH Standard Foot Shell, With CAP	FS	Sold separately (see catalogue)
RUSH Sandal Toe Foot Shell, No CAP	FTC-4R-1	
Black Spectra sock	SS	Included
RUSH EVAQ8 Rebuild Kit	EVRB	Sold separately
RUSH EVAQ8 Release Valve	EVRV	Sold separately
RUSH ROGUE 2 H2O		
Part description	Part number	Included / Sold separately
RUSH ROGUE 2 H2O Foot with regular toe	H2R2-xx-y-RU	Included

*xx = Size, y = Category (see table §3).

2. DESCRIPTION

The RUSH ROGUE 2 is an energy-return prosthetic foot with vertical compression and torsional rotation. It consists of:

- A fiberglass keel
- A fiberglass sole plate
- A rubber heel shock bumper
- A male pyramid connection

The RUSH ROGUE 2 EVAQ8 incorporates an elevated vacuum system inside the rubber heel shock bumper. The RUSH ROGUE 2 H2O incorporates an integrated Vibram® sole and custom molded rubber toe-wrap for superior traction on slick surfaces.

RUSH ROGUE 2 Foot and RUSH ROGUE 2 EVAQ8 are available in regular and right or left sandal toe versions and are delivered with a Spectra sock and a set of heel wedges.

3. PROPERTIES

Foot	RUSH ROGUE 2 Foot	RUSH ROGUE 2 EVAQ8	RUSH ROGUE 2 H2O*
Weight*	1039 g / 2.29 lbs	1048 g / 2.31 lbs	883 g / 1.96 lbs
Build height*	22 to 24 cm	155 mm / 6.13"	142 mm / 5.63"
	25 to 27 cm	168 mm / 6.63"	152 mm / 6.00"
	28 to 29 cm	171 mm / 6.75"	165 mm / 6.50"
Heel height	10 mm / 3/8"		

*Weight based on size 26, cat. 4, with foot shell and Spectra sock

Build height based on sizes 23, 26 or 29, cat. 4, with foot shell, Spectra sock, and 10 mm heel height

H2O weight and build height does not include foot shell or Spectra sock

This device has been tested according to ISO 10328 for a maximum patient weight up to 166 kg for 2 million cycles, corresponding to a service life of 2 to 3 years depending on patient's activity level.

Selection of foot category based on patient's weight and activity level										
Weight*)	lbs	0-105	106-140	141-175	176-210	211-245	246-280	281-315	316-350	351-365
	kg	0-48	49-64	65-79	80-95	96-111	112-127	128-143	144-159	160-166
Activity Level	Moderate	1	2	3	4	5	6	7	8	9
	High	1	2	3	4	5	6	7	8	9
	High Impact	2	3	4	5	6	7	8	9	-

*) Body mass limit not to be exceeded (ISO 10328)

4. MECHANISM OF ACTION


The half rocker sole plate provides constant contact to the ground while eliminating dead spots, providing a smooth rollover, and natural energy return. The heel shock bumper stores energy during initial contact and loading response which is then released during midstance. The instep absorbs energy during terminal stance and releases the stored energy during pre-swing.

5. INTENDED USE/INDICATIONS

This medical device is supplied to healthcare professionals (prosthetists) who will train the patient in its use. The prescription is made by a doctor who assess the patient's ability to use the device.

⚠ This device is for multiple use on a **SINGLE PATIENT**. It must not be used on another patient.

This device is intended to be integrated in a custom-made external lower limb prosthesis to ensure the function of the foot in patients with unilateral or bilateral lower limb amputation or deficiencies (transtibial/transfemoral amputation, knee/hip disarticulation, congenital limb deficiencies). It's intended for patients who would benefit from the dynamic toe response, vertical compression, and torsional rotation.

This device is indicated for patients with moderate to high activity level (K3 to K4).
Maximum weight (load carrying included): 166 kg / 365 lbs (See table §3)

6. CLINICAL BENEFITS

- Walking comfort
- Shock absorption and reduction of socket forces
- Axial rotation
- Possible ambulation on uneven ground

7. ACCESSORIES AND COMPATIBILITY

For RUSH ROGUE 2 FOOT and RUSH ROGUE 2 EVAQ8, an appropriate foot shell must be mounted on the foot (refer to our catalogue). The foot includes a male pyramid connection designed to be compatible with standard female pyramid connectors and adapters (refer to our catalogue).

8. MOUNTING

After dynamic alignment, tighten the pyramid adjustment screws according to the specifications of the connector manufacturer. Secure pyramid adjustment screws with a thread locking adhesive (i.e., Loctite 242).

RUSH ROGUE 2 EVAQ8

The straight barb, exhaust filter, tubing, inline filter, socket right angle barb and housing, tube securing hook and loop tape are included with the foot and will require assembly prior to use.

To connect the vacuum system to the socket:

- Connect a short piece of the tubing to the socket right angle barb.
- Install the inline filter into the distal end of the tubing (the inline filter may be placed anywhere in the tubing run between the socket and vacuum pump heel).
- Connect a second piece of tubing to the distal end of the inline filter and route the tube to the medial side of the pylon or wrap the tubing around the pylon (to prevent damage to the tube or snagging while walking).
- Cut the tubing to desired length and connect to the straight barb located in the recessed area of the vacuum pump heel.
- Secure the tubing to the pylon using the included hook and loop tape or other appropriate tape.

Spectra® sock

A Spectra sock is included to protect the foot shell and fiberglass components and minimize noise. It must be placed over the keel and sole plates before mounting the foot shell.

NOTE: The RUSH ROGUE 2 H2O is not intended to be worn with a foot shell. Therefore, it does not include a Spectra sock.

Foot shell

To install and remove the foot shell, use a shoehorn to prevent damage of the foot module.

⚠ Never remove the foot from the foot shell by pulling manually. Never use a screwdriver or any other inappropriate instrument to remove it. This could damage the foot.

⚠ The RUSH ROGUE 2 H2O is not intended to be worn with a foot shell. Doing so may damage the foot module.

9. ALIGNMENTS

Bench alignment

Recommendation: Add a 10 mm (3/8") wedge under back of the foot during bench alignment.


Sagittal plane:

- Define the appropriate socket flexion and heel height
- Position the load line so that it falls through the center of the foot pyramid connection. As the fiberglass composite material provides substantially more flexibility than the other prosthetic feet, this represents a good starting point for a bench

alignment

Frontal plane:

- Define the appropriate socket adduction/abduction
- Position the load line so that it falls through the midline of the foot in a neutral mediolateral position

Static alignment

Due to the shape of the rocker sole (bottom blade), patients may experience a new sensation when seeking the midpoint of the foot. The rocker sole allows patients to find their own comfortable static or standing position. The rocker bottom shape of the sole serves 2 primary functions:

- Provide a continuous and progressive point of contact throughout the entire step.
- Eliminate any “flat” or “dead” spot.

At this stage of alignment, adjustment of the anteroposterior set screws at the proximal connector is the more appropriate place to make changes regarding plantarflexion or dorsiflexion, rather than the introduction of a wedge.

Dynamic alignment

The use of a slide adapter is recommended for dynamic alignment purposes, as it provides the best solution to troubleshooting the most common alignment issues listed below:

- Hard or soft heel
- Hard or soft toe
- Varus or valgus movements during stance phase

Once the optimal relative socket/foot alignment is determined, use plantarflexion or dorsiflexion at the proximal foot adapter to optimize heel-to-toe comfort and energy return.

10. ADJUSTMENTS

If the patient still requests additional heel stiffness, use the heel wedges provided. These heel wedges use a double-sided sticker so the lower face of the sole plate should be degreased before using. See figures below for correct wedge placement.

NOTE: The RUSH ROGUE 2 H2O does not include heel wedges.

CORRECT HEEL WEDGE PLACEMENT	INCORRECT HEEL WEDGE PLACEMENT
	
<p>Install Heel Wedge on plantar surface of rocker approximately 3.2 mm (1/8") forward from posterior end.</p>	<p>Placement between components will adversely affect foot performance.</p>

11. TROUBLESHOOTING

- ⚠ If the patient notices any abnormal behavior or feel any changes in the characteristics of the device (noise, play, excessive wear...), or if the device has received a severe impact, he should stop using the device and consult its prosthetist.

12. WARNINGS

- ⚠ In case of damaged packaging, check the integrity of the device.
- ⚠ Never use the foot without a foot shell and a Spectra sock (except for RUSH ROGUE 2 H2O).
- ⚠ Never loosen the pyramid fastening screws.
- ⚠ The patient must inform its prosthetist if he gains or loses weight significantly.
- ⚠ Make sure that the foot and inside of the foot shell are free of impurities (e.g., sand). The presence of impurities causes the fiberglass parts and the foot shell to wear out. Clean the foot according to the instructions (see §17).
- ⚠ After swimming or using in water, the foot including its foot shell must be cleaned (see §17).
- ⚠ Failure to follow the instructions for use is dangerous and will void the warranty.

13. CONTRAINDICATIONS

- ⚠ Use for a patient whose maximum weight (load carrying included) may exceed 166 kg / 365 lbs
- ⚠ Use for a patient who do not meet the requirements of an activity level of K3 or higher
- ⚠ Use for activities associated with a risk of significant impact or excessive overloading

14. SIDE EFFECTS

There are no known side effects directly associated with the device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

15. MAINTENANCE AND CONTROL

No maintenance operation such as lubrication, work on the screws or other parts is required.




The foot module must be inspected by the prosthetist at least every six months. Inspections at shorter intervals are required if the user is more active.

The Spectra sock and the foot shell must be replaced by the prosthetist at regular intervals, depending on the patient’s level of activity. If these parts are damaged, it can lead to premature foot wear.

The **RUSH FOOT EVAQ8** components (tubing, inline filter, one-way valves housed inside the vacuum heel) may need periodic replacement during the life cycle of the system and are not replaceable under the warranty as it is considered normal wear.

16. PERIODIC INSPECTION OF THE RUSH ROGUE 2 EVAQ8 SYSTEM

- Visually inspect the tubing for kinks, cracks, or wear that may leak air into the system. Replace tubing if any of these conditions exist.
- Remove the inline filter from tubing and look through it. If light can be seen, the filter is clean. If light is blocked, blast air from a syringe through the inline filter from distal to proximal end (reverse of normal flow) to attempt to clear the blockage. If blockage persists, the filter needs to be replaced.
- The one-way valves contained in the vacuum heel may need to be cleaned and flushed with distilled water or isopropyl alcohol to ensure proper function. This procedure should be done only by a qualified professional.
- To flush the one-way valves and vacuum heel:

					
1. Exhaust Assembly 2. Exhaust Filter 3. Large O-Ring 4. Valve Body Adapter 5. Small O-Ring	2. Exhaust Filter 6. Duckbill Valve	7. Straight Barb	8. Duckbill Valve	9. Paperclip (foot on side)	9. Paperclip (foot upright)

- Remove the vacuum hose from the user’s socket keeping it attached to the EVAQ8 foot.
- Remove the EVAQ8 foot from the user’s socket.
- Remove the vacuum hose from the EVAQ8 foot.
- Using a 5/16” socket, remove the (1) Exhaust Assembly from the heel, the (4) Valve Body Adapter will most likely stay attached to the (2) Exhaust Filter.
NOTE: If (4) Valve Body Adapter does not remove with (2) Exhaust Filter, use a flat blade screwdriver to remove it.
- Place the (4) Valve Body Adapter into a soft jaw vice or grip with soft face pliers and remove the (2) Exhaust Filter with 5/16” socket.
NOTE: (6) Duckbill Valve will be slightly pressed into the bottom of the exhaust filter.
- Using a 1/4” deep socket, remove the (7) Straight Barb from the other side of the heel.
- Inside the heel beneath where the (7) Straight Barb was located is another (8) Duckbill Valve. Remove the (8) Duckbill Valve by tapping the foot against your hand or by straightening out a (9) Paperclip and inserting it into the other side of the heel to push out the (8) Duckbill Valve.
- Inspect both (3 & 5) O-rings on the (4) Valve Body Adapter. There is one at the base of the threads, and one in a groove on the body. Replace both if they show any wear.
- Clean the female threads in both sides of the heel with a cotton swab and isopropyl alcohol or distilled water.
- If you are reusing the (1) Exhaust Filter, (7) Straight Barb, in-line filter, and (8) Duckbill Valve, clean them with isopropyl alcohol or distilled water. Take extra care and make sure the (8) Duckbill Valve are clean and free from debris (a magnifier is helpful for inspection). Flush the in-line filter from both directions to ensure it is clean.
- Once the parts are dry, or if you are using new parts from a rebuild kit, set all the parts out on a clean surface.
- Insert the turquoise (8) Duckbill Valve into the threaded opening of the (1) Exhaust Filter so that the (8) Duckbill Valve rim is flush against the opening of the (1) Exhaust Filter and the tip of the (8) Duckbill Valve is inside the (1) Exhaust Filter. Thread the (1) Exhaust Filter into the (4) Valve Body Adapter **BY HAND** until snug.
- **BY HAND**, thread the assembled (1) Exhaust Assembly into side of the heel that will **NOT** be connected to the prosthetic

socket by vacuum hose.

- Once the (1) Exhaust Assembly is threaded in snugly **BY HAND**, torque to 15 in-lbs. Do not over torque. Over torquing will snap the threads and will not be covered under warranty.
NOTE: If you do not have a torque wrench, thread the (1) Exhaust Assembly until you feel a hard stop and then turn 1/16 of a revolution more.
- Insert the white (8) Duckbill Valve into the recessed side of the heel so the tip of the (8) Duckbill Valve points into the heel. Using a small screwdriver or a straightened-out (9) paperclip, make sure the (8) Duckbill Valve is seated all the way into the recess.
- BY HAND**, thread the (7) Straight Barb into the side of the heel that is recessed.
- Once the (7) Straight Barb is threaded in snugly **BY HAND**, torque it to 15 in-lbf. This is a very low torque value, and over torquing will snap the threads on the (7) Straight Barb and will not be covered under warranty.
- If you do not have a torque wrench, thread the (7) Straight Barb until you feel a hard stop and then turn 1/16 of a revolution more.
- Reattach vacuum hose to EVAQ8 foot by sliding the vacuum hose over the (7) Straight Barb.
- Put the Spectra® Sock and foot shell over the EVAQ8 foot.
- Reattach the EVAQ8 foot to the user's socket.
- Reattach the other end of the vacuum hose to the user's socket. The vacuum hose may be routed to the prosthetist's preference.

17. CLEANING

Remove the foot shell and the Spectra sock, rinse the foot in clear water, clean it with neutral soap and dry carefully. The foot shell can be cleaned with a damp cloth or sponge. It must be dry before reuse.

The device is not resistant to solvents. Exposure to solvents may cause damage.

18. ENVIRONMENTAL CONDITIONS

Temperatures range for use and storage: -20 to 60°C [-4 to 140°F]

Relative air humidity: no restrictions

Waterproof: the device is resistant to fresh, sea, and chlorinated water.

The foot shell is not resistant to ultraviolet light (UV). Do not store it under direct sunlight.

19. DISPOSAL

The foot is made of fiberglass composite material, rubber, Vibram (H2O), metal and plastic (EVAQ8). The foot shell is made of thermosetting polymer material. The device and its packaging must be disposed of in accordance with local or national environmental regulations.

20. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of 1 st declaration
	Authorized representative in the European Union		Single patient, multiple use		

21. REGULATORY INFORMATION

This product is a CE-marked medical device and is certified as conforming with Regulation (EU) 2017/745.

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PROTEOR SAS