RUSH Chopart®

CHO

- | | | | |
|-----------|--------------------------|-----------|-----------------------------|
| FR | Notice d'Utilisation | FI | Käyttöohjeet |
| EN | Instruction for use | PL | Instrukcja użytkowania |
| DE | Gebrauchsanweisung | CS | Návod k použití |
| IT | Istruzioni per l'uso | SK | Návod na používanie |
| ES | Instrucciones de uso | HR | Upute za uporabu |
| PT | Instruções de utilização | RU | Инструкция по использованию |
| NL | Gebruikershandleiding | UK | Інструкція з використання |
| DA | Brugervejledning | JA | 取扱説明書 |
| NO | Bruerveiledning | ZH | 使用说明 |
| SV | Bruksanvisning | AR | تعليمات الاستخدام |

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RUSH Chopart®
Instructions for use for prosthetists
Read before use

IFU-02-014
 Rev. A
 2023-08

Pass on § 10 to 17 of these instructions to the patient.

1. INCLUDED ITEMS

Part description	Part number	Included / Sold separately
RUSH Chopart	CHO-xx-x-RU*	Included
Spectra sock	S0-NPS-200xx-00*	Suitable sock included
Foot shell	FS-xx-xx*	Sold separately

* see catalog

2. DESCRIPTION

RUSH Chopart is an energy-return prosthetic foot plate composed of fiberglass composite used for the manufacture of partial-foot prostheses.

3. PROPERTIES

Weight*	380 g / 0.84 lbs		
Size	22 to 24 cm	25 to 27 cm	28 to 30 cm
Build height*	25 mm / 1.00"	29 mm / 1.13"	22 mm / 0.875"
Heel height	10 mm / ³ / ₈ "		

*Weight based on size 26, cat. 4, with foot shell and Spectra sock

Build height based on sizes 23, 26 or 29, cat. 4, with foot shell, Spectra sock, and 10 mm heel height

Selection of foot category based on patient's weight and impact level

Weight	lbs	0-105	106-140	141-175	176-210	211-245	246-280	281-315	316-350	351-365
	kg	0-48	49-64	65-79	80-95	96-111	112-127	128-143	144-159	160-166
Impact level	Low	1	2	3	4	5	6	7	8	9
	Moderate	1	2	3	4	5	6	7	8	9
	High	2	3	4	5	6	7	8	9	-

4. MECHANISM OF ACTION

The curved sole plate provides constant contact to the ground while eliminating dead spots, providing a smooth rollover, and natural energy return.

5. INTENDED USE/INDICATIONS

This medical device is supplied to healthcare professionals (prosthetists), who will train the patient in its use. The prescription is drawn up by a doctor who assesses the patient's ability to use the device.

⚠ This device is for multiple use on a **SINGLE PATIENT**. It must not be used on another patient.

This device is intended to be integrated in a custom-made partial-foot prosthesis to ensure the function of the foot in patients with unilateral or bilateral partial-foot amputation or deficiencies (Chopart, Lisfranc, Pirogoff or Boyd amputations, congenital foot deficiencies). It's intended for patients who would benefit from the dynamic toe response.

This device is indicated for patients with medium to high activity level (K2 to K4) for walking and physical activities without excessive overload.

Maximum weight (load carrying included): 166 kg / 365 lbs (See table §3)

6. CLINICAL BENEFITS

- Possible ambulation on uneven ground
- Smooth roll-over


7. ACCESSORIES AND COMPATIBILITY

The RUSH Chopart must be bonded and laminated onto the distal end of a prosthetic socket for partial foot amputation. An appropriate foot shell must be mounted on the foot (refer to our catalog).

8. ALIGNMENT

Recommendation: Add a 10 mm ($\frac{3}{8}$ ") wedge under back of the foot during bench alignment.

- Integrate the foot plate to the socket with an initial bench alignment.
- Adjust the heel height appropriate to the shoe.
- Align the plate on the socket by setting the plantar flexion / dorsiflexion and inversion / eversion.
- The load line of the socket should align with the foot plate at 2/3 of the posterior distance of the plate.
- Set the rotation of the foot plate.

9. ASSEMBLY

Socket initial integration

Integrate the bench aligned RUSH Chopart plate into the test socket by following the steps below. It is recommended to bond only the posterior 1/3 of the foot plate to the socket.

- Clean upper and lower surfaces of the foot plate thoroughly with clean cloth and acetone to remove all release agents, oils, etc.
- Follow proper air ventilation and skin protection as required.
- Lightly scuff the upper and lower bonding surfaces with fine grit paper just enough to remove the sheen off the surface.
- Avoid deep sanding of the surface to protect the fibers of the laminate.
- Epoxy the foot plate to the socket by filling in the edges and providing a smooth transition.
- Optimize the dynamic performance of the plate by adhering as closely as possible to potting 1/3 of the plate.
- Secure the foot plate with synthetic casting tape, scotch cast or other temporary wrap materials.
- Perform dynamic alignment and adjust the plate as required.

Socket lamination

- Secure the orientation of the dynamic alignment procedure using an alignment device or jig.
- Prepare for direct lamination by removing temporary bonding materials and synthetic casting tape.
- Prepare the RUSH Chopart plate for final lamination by cleaning the exposed surfaces and ensuring surface preparations with a light scuff using fine grit paper.
- Laminate the plate to the socket with carbon fiber cloth starting with a long posterior layer and follow with medial / lateral layers.
- Fold the long posterior layer back over medial / lateral layers to interlock the lamination sequence.
- Secure with circumferential wraps.
- Compensate for the weight and activity level of the user with increased laminating layers for additional strength.

Spectra sock

A Spectra sock is included to protect the foot shell and fiberglass components and minimize noise. It must be placed over the foot before mounting the foot shell.

Foot shell

To install and remove the foot shell, use a shoehorn to prevent damage of the foot module.

- ⚠ Never remove the foot from the foot shell by pulling manually. Never use a screwdriver or any other inappropriate instrument to remove it. This could damage the foot.

10. TROUBLESHOOTING

- ⚠ If the patient notices any abnormal behavior or feel any changes in the characteristics of the device (noise, play, excessive wear...), or if the device has received a severe impact, they should stop using the device and consult their prosthetist.

11. WARNINGS

- ⚠ In case of damaged packaging, check the integrity of the device.

- ⚠ Never use the foot without a foot shell and a Spectra sock.
- ⚠ The patient must inform their prosthetist if they gain or lose significant weight.
- ⚠ Make sure that the foot and inside of the foot shell are free of impurities (e.g., sand). The presence of impurities causes the fiberglass parts and the foot shell to wear out. Clean the foot according to the instructions (see §15).
- ⚠ After swimming or using in water, the foot including its foot shell must be cleaned (see §15).
- ⚠ Failure to follow the instructions for use is dangerous and will void the warranty.

12. CONTRAINDICATIONS

- ⚠ Do not use for a patient whose maximum weight (load carrying included) may exceed 166 kg / 365 lbs.
- ⚠ Do not use for a patient who do not meet the requirements of an activity level of K2 or higher.
- ⚠ Do not use for activities associated with a risk of significant impact or excessive overloading.

13. SIDE EFFECTS

There are no known side effects directly associated with the device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

14. MAINTENANCE AND CONTROL

No maintenance operation such as lubrication, work on the screws or other parts is required.

The foot module must be inspected by the prosthetist at least every six months. Inspections at shorter intervals are required if the user is more active.

The Spectra sock and the foot shell must be replaced by the prosthetist at regular intervals, depending on the patient's level of activity. If these parts are damaged, it can lead to premature foot wear.

15. CLEANING

Remove the foot shell and the Spectra sock, rinse the foot in clear water, clean it with neutral soap and dry carefully. The foot shell can be cleaned with a damp cloth or sponge. It must be dry before reuse.

- ⚠ The device is not resistant to solvents. Exposure to solvents may cause damage.

16. ENVIRONMENTAL CONDITIONS

Temperature range for use and storage: -20°C to 60°C [-4 to 140° F]

Relative air humidity: no restrictions

Waterproof: the devices are resistant to fresh, sea, and chlorinated water.

- ⚠ The foot shell is not resistant to ultraviolet light (UV). Do not store it under direct sunlight.

17. DISPOSAL

The device is made of fiberglass composite material. The device and its packaging must be disposed of in accordance with local or national environmental regulations.

18. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of 1 st declaration
	Authorized representative in the European Union		Single patient, multiple use		

19. REGULATORY INFORMATION

This product is a CE-marked medical device and is certified as conforming with Regulation (EU) 2017/745.

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PROTEOR SAS