
## KINEGEN.wave


Polyzentrisches Kniegelenk  
Polycentric Knee Joint  
Rodilla policéntrica

3A900      3A900/1


1




2


3



7

8

9

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## Content

1	Technical Data.....	20
2	Description and Function.....	20
2.1	Application .....	20
2.2	Construction.....	20
3	Assembling .....	20 – 22
3.1	Basic set-up.....	21
3.2	Positioning of the socket.....	21
3.3	Limitation of flexion .....	22
3.4	Replacing the extension Spring.....	22
4	Adjustment .....	22 – 24
4.1	Adjustment of the standing safety.....	22
4.2	Adjustment of the swing phase.....	23
4.3	Individual adjustment of the lock.....	24
4.4	Individual adjustment of the extension stop (Fig. 4) .....	24
5	Finishing the prosthesis.....	25
6	Technical Notes .....	25
7	Maintenance and service .....	25
8	Spare parts.....	25
9	Use .....	26
9.1	Restriction of Use.....	26
9.2	Misuse of Product .....	26
9.3	Useful Life.....	26
10	Guarantee.....	27 – 28
10.1	Conditions .....	27
10.2	Guarantee and Preclusion period .....	27
10.3	Scope of Guarantee and Exclusion .....	28
10.4	Guarantee service .....	28
11	Declaration of Conformity.....	29
12	Storage Conditions .....	29
13	Disposal .....	29

# 1 Technical Data

	3A900	3A900/1
Connection proximal:	Adjustment Core	M36 thread

The **3A900** and **3A900/1 KINEGEN.wave** polycentric knee joints are approved for **activity level 1 and 2** according to Streifeneder walky-system and for maximum body weight of up to **125 kg / 275 lbs**, including carried and supported loads.

## 2 Description and Function

### 2.1 Application

The products are to be used exclusively for lower limb prosthetics.

### 2.2 Construction

The 3A900 and 3A900/1 KINEGEN.wave knee joints are polycentric knee joints with an integrated extension assist and an optional lock that can be released by a release cord. The upper and lower parts of the four-axis knee joint are connected by levers. The integrated extension assist is used to control the swing phase. The knee joint is connected proximally by means of an Adjustment Core (3A900) or M36 thread (3A900/1), distally by means of a tube holder with tube clamp (diameter)  $\varnothing$  30 mm.

## 3 Assembling

The arrangement of the knee joint axes ensures that in the extended position the momentum pivot point (D) is above the knee joint and behind the load line (Fig. 1). Therefore, the stance phase safety is determined exclusively by the prosthesis alignment.

### Attention!

Only with a correct prosthetic alignment can the properties of the knee joint be used optimally.

For correct assembly, the residual limb position must be taken into account for positioning the socket connection.

To facilitate final positioning of the socket adapter, plumb lines can be drawn in the frontal plane and sagittal plane from the hip joint pivot point during cast removal and fitting of the test socket (see Fig. 2 and Fig. 3).

### **Attention!**

When bending the knee joint, do not reach into the joint mechanism – risk of injury by pinching the fingers!

## **3.1 Basic set-up**

The prosthesis is set up in such a way that in frontal view the plumb line (a) runs from the centre of the socket over the Adjustment Core and the centre of the tube holder downwards. In sagittal view, the construction reference line (b) runs from the greater trochanter (T) over the anterior upper joint axis to the transition between the posterior third and middle third of the prosthetic foot, the heel height is set at approx. 5 mm (Fig. 2).

## **3.2 Positioning of the socket**

Make a temporary mark (T and M) on the lateral side proximally and distally in the centre of the prosthetic socket and connect the mark with a line (c) from the upper edge of the socket to the end of the socket. The connecting line (c) must intersect with the construction reference line in the area of the socket entry plane (Fig. 3). The basic setting of the socket flexion is normally 5°, whereby the individual situation (for example hip joint contracture) must be taken into account.

### **Attention!**

If residual limb flexion is not taken into account, the prosthesis wearer may not be able to actively extend the knee joint. This can lead to malfunctions or reduced stability.

Shifting the shaft posteriorly increases the dynamics and reduces the stability. Shifting the shaft anteriorly increases the stability and reduces the dynamics. The connection between the prosthetic socket and the knee joint must be made in such a way that the knee joint is rotated outwards by about 5°. Abduction and adduction can be performed using the adjustment screws on the socket adapter and prosthetic foot. The fine adjustment takes place during the dynamic fitting.

### 3.3 Limitation of flexion

The KINEGEN.wave enables a knee flexion of 145°. However, this flexion angle can only be achieved with a correct socket/prosthesis alignment.

Clamps of socket adapters must be located medially or laterally on the shaft according to the instructions for use. If the clamp is posterior, the flexion angle of the knee joint may not be fully utilised or the knee joint may be damaged.

With a correct prosthetic alignment, it can be assumed that the user's posterior touches the prosthetic foot before the socket comes into contact with the knee joint.

### 3.4 Replacing the extension Spring

To replace the spring, turn screw C (Fig. 9) anticlockwise until it can be removed from the knee joint. Replace the old spring from the knee joint with the new 3A900/E1 spring. Before fitting, make sure that the spring has been greased.

#### Attention!

Note that when using the replacement spring, the default setting of screw c is as shown in Figure 9, with the screw protruding by 3 mm. If necessary, the screw can be turned in a maximum of one revolution to achieve a faster extension in the swing phase.

## 4 Adjustment

### 4.1 Adjustment of the standing safety

#### Attention!

If the prosthetic socket is tilted more anteriorly at the proximal pyramidal connection of the knee joint, standing safety is increased. If the prosthetic socket tilts posteriorly, the dynamics increase and the standing safety decreases.

#### Attention!

Standing safety is also influenced by the construction of the prosthetic foot.

## 4.2 Adjustment of the swing phase

The combination of axial friction and the extension spring of the knee joint is used to adjust the extension and flexion movement of the knee during the swing phase. When delivered, the KINEGEN.wave is set so that the extension movement of the knee joint is not impeded. The extension assist supports the extension of the knee joint and counteracts flexion. Optimal adjustment of the swing phase enables the prosthesis wearer to have a harmonious gait pattern and an energy-efficient gait. The swing phase behaviour and thus the gait pattern is influenced, among other things, by the lower leg length, the weight of the foot/shoe and the activity of the user.

The sequence for adjusting the swing phase when delivered is the following (Fig. 4).

1. By turning screw A clockwise, the axle friction and thus the resistance to movement is increased. To reduce the heel lift again, turn the adjustment screw A anticlockwise.
2. By turning adjustment screw C, the force of the extension assist is adjusted to synchronise the swing motion of the prosthesis to that of the healthy side. To allow a faster forward swinging motion, turn adjustment screw C clockwise. After each readjustment or adjustment, the axle friction must be readjusted as described in point 1.

### **Attention!**

If the strength of the standard spring for extension assist is not sufficient, it can be replaced with a 3A900/E1 spring, as described in 3.4.

### 4.3 Individual adjustment of the lock

#### **Adjustment of the Release Cord (Fig. 5 and Fig. 6)**

The Release Cord must be shortened individually for the patient. To do this, loosen the two set screws on the handle, adjust the length of the locking cable and retighten the set screws hand-tight. The locking cable must be guided through both channels attached to the handle (Fig. 5 C). Attach the handle to the outside of the socket. Make sure that the locking cable is not under tension when the joint is locked. At the same time, it must be ensured that the available pulling distance is sufficient to trigger the lock. If the release cable is correctly set, the user can unlock the knee joint for sitting down via the manual operation (Fig. 6).

#### **Temporarily deactivating the lock (Fig. 7)**

The user is able to walk with a polycentric knee joint at times, e.g. under supervision, otherwise he needs a locked knee joint. Remove one of the two lugs on the handle of the release cable to allow the handle to be hooked in when unlocked (Fig. 7). Make sure that the locking bolt cannot engage when the unit is hooked in. Thus, the knee can be used temporarily as a polycentric knee joint.

#### **Temporarily deactivating the lock (Fig. 8)**

The user is able to walk with a polycentric knee joint at times, e.g. under supervision, otherwise he needs a locked knee joint. Hold the lock open with the locking cable and screw in the securing screw E.

### 4.4 Individual adjustment of the extension stop (Fig. 4)

To adjust the extension stop individually during fitting, turn screw B with a hexagonal screwdriver. To ensure an earlier extension stop, the screw should be turned in clockwise.

After adjusting the extension damper, check that the lock continues to engage independently in full extension.

## 5 Finishing the prosthesis

### Attention!

The tightening torque for the adjusting screws is **14 Nm** with the 173P11 Torque Wrench, calibrated.

## 6 Technical Notes

### Attention!

Please take special care, that the 3A900 and 3A900/1 Kneejoint are not exposed to chemicals and substances such as acids, dust or sand!

**After each use in salt- or chlorinated (pool) water, the prosthesis must be thoroughly rinsed with fresh water. In addition, the prosthesis should be dried after rinsing.**

Failure to observe these instructions will make any claims against Streifeneder ortho.production GmbH null and void.

**Please inform your prosthesis wearer accordingly!**

## 7 Maintenance and service

### Attention!

The 3A900 and 3A900/1 knee joint must not be disassembled.

The 3A900 KINEGEN.wave knee joint must be checked for contamination and function after 6 months. The knee joint must be checked annually for wear and functionality and readjusted if necessary.

If necessary, clean the knee joint (not with compressed air) and adjust it again. Increased noise occur, please send the knee joint to Streifeneder ortho.production GmbH for maintenance. Excessive wear of the bump stops can impair their function.

## 8 Spare parts

Description	Item-No.
Alternative spring for extension assist	3A900/E1

## 9 Use

This product requires qualified and individual adaption to the patient according to medical indication. All adjustments must be done by qualified professional staff.

### 9.1 Restriction of Use

Because of the fitting and for hygienic reasons, we have to point out that the product is to be used only by one patient and must not be given to any other patient.

#### **Attention!**

We are pointing out, that this product is tested and released only in combination with all other Streifeneder modular parts.

### 9.2 Misuse of Product

This product is exclusively constructed for the activity level (walky) 1 to 2 up to maximum body weight, including carried or supported loads, of up to 125 kg and is not suitable for sporting activities (e.g. jumping, running, etc.).

#### **Attention!**

The approved maximum weight load of this product is limited and depends on the component with the lowest weight limit used in the prosthetic construction, which determines the maximum weight load of the entire construction. The maximum weight load means the body weight of the user including carried and supported loads is designed for up to 125 kg.

### 9.3 Useful Life

The product has been tested to ISO 10328 for three million load cycles. Depending on the activity of the prosthesis wearer, this corresponds to a duration of use of 36 to 60 months.

## **10 Guarantee**

The following guarantee bond shall apply in conjunction with any legal statutory warranty claims and completely independent thereof. It is provided in lieu of any other expressed or implied warranties or suitabilities for a particular purpose (to the extent permitted by law).

### **10.1 Conditions**

Claims resulting from this manufacturers guarantee can only be asserted, if the relevant product was installed or adjusted by or under the direction of a certified / licensed orthopedic technician AND if the product has been properly and regularly maintained by a trained technician within the provided intervals - as far as such maintenance and inspection intervals are predefined by Streifeneder ortho.production GmbH. The conducted maintenance and service checks must be confirmed with the date, official stamp and signature of the performing specialist company. The corresponding verifications shall be submitted to Streifeneder ortho.production GmbH without request, upon the assertion of claims under this Guarantee.

### **10.2 Guarantee and Preclusion period**

The duration of our manufacturers guarantee varies depending on the product. However, the duration of the Guarantee starts in any case upon delivery to the customer. The special terms and conditions of the Guarantee for this particular product are:

Guarantee Duration: 3A900 and 3A900/1 KINEGEN.wave polycentric Knee joints  
= 24 months

Claims under this Guarantee must be asserted immediately in writing to Streifeneder ortho.production GmbH, at the latest within one month after knowledge of the existence of the preconditions for claims under this Guarantee. Rights to assert guarantee claims are expressly excluded after the expiry of this time limit.

### **10.3 Scope of Guarantee and Exclusion**

It is guaranteed that the product is free from defective material and / or workmanship.

The Guarantee does not cover damages resulting from

- faulty installation,
- operating errors,
- inappropriate application (especially contrary to existing manufacturer specifications / instructions)
- normal wear and tear
- external influences such as force majeure or interventions by third parties.

Likewise not included are merely optical impairments without influence on the functionality of the relevant product, as well as defects arising from the installation of damaged components, the application of spare parts not approved by Streifeneder ortho.production GmbH and / or interventions or repair works conducted by distributors / specialist companies not authorized by Streifeneder ortho.production GmbH.

This Guarantee does not create a liability of Streifeneder ortho.production GmbH with regard to consequential loss or damage of any kind or other claims for compensation with the exception of mandatory legal provisions according to which Streifeneder ortho.production GmbH can be held accountable as manufacturer beyond the scope of this voluntary Guarantee.

### **10.4 Guarantee service**

If all preconditions of this Guarantee are met, Streifeneder ortho.production GmbH will, at its sole discretion, either repair or replace the product free of charge.

## 11 Declaration of Conformity

Streifeneder ortho.production GmbH declares as manufacturer in sole responsibility, that the 3A900 and 3A900/1 KINEGEN.wave polycentric Knee joints complies with the Regulation (EU) 2017/745 for medical devices (MDR).

## 12 Storage Conditions

Store the product in its original packaging at room temperature in a dry place.

## 13 Disposal

The product is to be recycled or disposed of without risk to humans or the environment. Please observe the symbol on the product label (packaging).

This instruction manual is protected by copyright.

Reproduction of the complete instruction manual, even in abstracts, in any form or by any means, is only permitted with explicit written authorisation of the Streifeneder ortho.production GmbH.

The indicated measurements and weights are guide values without obligation. During the delivery period, changes in design or form, deviations in colour as well as modifications of the score of delivery shall remain reserved on part of the manufacturer, as long as the purchased product is not substantially altered and as long as the changes are reasonable for the buyer.

Products are subject to change without prior notice after press date of this instruction manual, 2021-09-30.

Our general terms and conditions apply; they can be viewed anytime at [www.streifeneder.com/op](http://www.streifeneder.com/op).

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