




Aluminium 30mm Structural Components

| Part Number | Manufacturer | Description | Weight Limit |
|---------------|--------------|--|--------------|
| 1* NAAA210 | TL | 4-Hole Pyramid Adapter | 160kg |
| 2* NAAA211 | TL | 4-Hole Female Adapter | 160kg |
| 3* NAAA212 | TL | 4-Hole Rotatable Adapter - Note: This part number is for the 4 hole base plate only. It will take either NAAA214-01 or NAAA216-01, these must be ordered separately. | 160kg |
| 4* NAAA214-01 | TL | Rotatable Female Insert | 160kg |
| 5* NAAA216-01 | TL | Rotatable Male Pyramid Insert - Note: This part number is for the insert only. It can be used in either the 4 bolt adapter NAAA212 or the 3-Prong laminating adapter NAAASS237, these must be ordered separately. | 160kg |
| 6* NAAA209 | TL | 4-Hole Single Offset Pyramid Adapter - Offset by 5mm in one direction only. | 160kg |
| 7* NAAA208 | TL | 4-Hole Double Offset Pyramid Adapter - Offset by 5mm in two directions. | 160kg |
| 8* NAAA100 | TL | Adjustable Height Clamp Adapter - This clamp allows up to 25mm adjustment in height, with the two bolt clamp being sufficient to hold the position, unless there is likely to be any significant impact in the patient's gait, such as jumping or running. | 160kg |
| 9* NAAA135 | TL | Pylon with Integral Receiver - This 30mm pylon with integral female receiver can be used with either the NAAA100 or NAAA225 tube clamps, but at 420mm long, is suitable for most transtibial applications and can be cut down to suit. | 160kg |

*Note: The weight limit is 123kg for the NAAA aluminium components above for patients of K4 activity level.

Manufacturer: TL = Trulife



Aluminium 30mm Structural Components

| Part Number | Manufacturer | Description | Weight Limit |
|-------------|--------------|---|--------------|
| 1* NAAA221 | TL | 35mm Double-Ended Adapter | 160kg |
| 2* NAAA222 | TL | 45mm Double-Ended Adapter | 160kg |
| 3* NAAA223 | TL | 60mm Double-Ended Adapter | 160kg |
| 4* NAAA224 | TL | 75mm Double-Ended Adapter | 160kg |
| 5* NAAA200 | TL | 420mm Bonded Pylon | 160kg |
| 6* NAAA202 | TL | 420mm Aluminium Pylon Bonded to Male Adapter | 160kg |
| 7* NAAA225 | TL | 30mm Clamp Adapter | 160kg |
| 8* NAAA285 | TL | Foot Pyramid - This foot adapter is specifically for use with the Seattle Lightfoot, Natural or SACH, though these now come with an adapter already fitted. | 136kg |

Stainless Steel 30mm Structural Components

| Part Number | Manufacturer | Description | Weight Limit |
|--------------|--------------|--|--------------|
| 1* NAAASS237 | TL | 3-Prong Laminating Piece Stainless Steel - This part number is only for the adapter and can take either NAAA214-01 or NAAA216-01, or any component with an M36 threaded top. | 136kg |

*Note: The weight limit is 123kg for the NAAA aluminium components above for patients of K4 activity level.

Manufacturer: TL = Trulife