






**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Titanium College Park Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	N4HPR-T	CP	4-Hole Pyramid Receiver, Standard	136kg
2	N4HPRR-T	CP	4-Hole Pyramid Receiver, Rotatable	136kg
3	N4HPA-T	CP	4-Hole Pyramid Adapter, Standard	136kg
4	N4HPAR-T	CP	4-Hole Pyramid Adapter, Rotatable	136kg
5	N4HPA-AO-T	CP	4-Hole Pyramid Offset, Axial (Offset by 7mm)	136kg
6	N4HPA-DO-T	CP	4-Hole Pyramid Offset, Diagonal (Offset by 7mm)	136kg
7	N4HPH-T	CP	4-Hole Pyramid Adapter with Pin-Lock Hole	136kg
8	N4HPHR-T	CP	4-Hole Pyramid Adapter with Pin-Lock Hole, Rotatable	136kg
9	NTA-T250-30	CP	Low-Profile Tube Adapter, 250mm, 30mm diameter	136kg
10	NTA-T430-30	CP	Low-Profile Tube Adapter, 430mm, 30mm diameter	136kg
11	NTC-T-30	CP	Tube Clamp, 30mm diameter	136kg
12	NDR-T35	CP	Double-Ended Adapter, 35mm	136kg

Manufacturer: CP = College Park



**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Titanium College Park Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
13	NDR-T45	CP	Double-Ended Adapter, 45mm	136kg
14	NDR-T60	CP	Double-Ended Adapter, 60mm	136kg
15	NDR-T75	CP	Double-Ended Adapter, 75mm	136kg
16	NTI-TTK-PA	CP	Transtibial Kit with Pyramid Adapter (30mm Tube Adapter - 250mm length, 30mm Tube Clamp, 4-Hole Pyramid Adapter)	136kg
17	NTI-TTK-PAH	CP	Transtibial Kit with Pyramid Adapter and Pin-Lock Hole (30mm Tube Adapter - 250mm length, 30mm Tube Clamp, 4-Hole Pyramid Adapter with Pin-Lock Hole)	136kg
18	NTI-TTK-PAR	CP	Transtibial Kit with Rotatable Pyramid Adapter (30mm Tube Adapter - 250mm length, 30mm Tube Clamp, 4-Hole Rotatable Pyramid Adapter)	136kg
19	NTI-TTK-PHR	CP	Transtibial Kit with Rotatable Pyramid Adapter and Pin-Lock Hole (30mm Tube Adapter - 250mm length, 30mm Tube Clamp, 4-Hole Rotatable Pyramid Adapter with Pin-Lock Hole)	136kg
20	NTI-TFK-PA	CP	Transfemoral Kit with Pyramid Adapter (30mm Tube Adapter - 430mm length, 30mm Tube Clamp, 4-Hole Pyramid Adapter)	136kg
21	NTI-TFK-PAH	CP	Transfemoral Kit with Pyramid Adapter and Pin-Lock Hole (30mm Tube Adapter - 430mm length, 30mm Tube Clamp, 4-Hole Pyramid Adapter with Pin-Lock Hole)	136kg
22	NTI-TFK-PAR	CP	Transfemoral Kit with Rotatable Pyramid Adapter (30mm Tube Adapter - 430mm length, 30mm Tube Clamp, 4-Hole Rotatable Pyramid Adapter)	136kg
23	NTI-TFK-PHR	CP	Transfemoral Kit with Rotatable Pyramid Adapter and Pin-Lock Hole (30mm Tube Adapter - 430mm length, 30mm Tube Clamp, 4-Hole Rotatable Pyramid Adapter with Pin-Lock Hole)	136kg

Manufacturer: CP = College Park