

RevoFit® Lamination Kit

The RevoFit Lamination Kit allows you to create built-in, user-controlled adjustability anywhere on a socket, allowing users make incremental changes to their socket fit and volume throughout the day without the need to add or remove socks.

For use with all commonly used lamination materials, the RevoFit uses a tried and tested system of manufacture. This kit has been formulated to resist the repetitive stress found in an adjustable socket, incorporating a high durability (HD), fatigue-resistant lace. The system is well established and empowers users by allowing adjustability; it is complemented by the user-friendly Click Reel providing additional benefits including a balance of precise adjustment and power, offering up to 7:1 mechanical advantage.

Waterproof*

*Waterproof note: Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

- Adjustable global and/or targeted compression
- Micro/macro adjustments to increase and to decrease socket volume
- Easy adjustments through clothing
- Optional grip rings to suit user dexterity needs
- Rapid donning and doffing for challenging residual limb shapes
- User Weight Limit: 135kg



Spend less time on remakes and refits.



Build adjustability, anywhere.

Structural and Socket Components Click Medical - RevoFit® Lamination Kit

Ordering Information

| RevoFit® Lamination Kit Part No. | Includes |
|----------------------------------|---|
| NCMRF-200-07-01 | <ul style="list-style-type: none"> • 1x Lamination Tube (1.8m) • 1x Reel Tool • 1x Lam. Collar |

Note: Please see page 308 for RevoFit® Lamination Kit spare parts ordering information.

Technical Information

| | | | |
|-----------------------------|----------------------------------|-------------------|------------------------------|
| RevoSurface® Tool Kit | | Lace | 2m HD (High Durability) Lace |
| Build Height of Click® Reel | 14mm visible profile, 24mm total | User Weight Limit | 300lb/135kg |
| Reel Diameter | 42mm | | |
| Reel Weight | 71g | | |
| Collar Diameter | 42mm | | |
| Tube/Sheath Length | 1.8m | | |
| Tube Outer Diameter | 3.2mm | | |

