



Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Special Adapters				
Part Number	Manufacturer	Description	Weight Limit	
1	NMA11A15/A	SF	Aluminium Tube Clamp Adapter with Pyramid Ø 30mm	100kg
2	NMA10A41/T	SF	Titanium Female Lamination Adapter with Shift	100kg
3	NMA10A43/T60	SF	60mm Titanium Double Adapter with Shift	100kg
	NMA10A43/T75		75mm Titanium Double Adapter with Shift	100kg
4	NMA11A14/A	SF	Aluminium Double Adapter with Pyramid	100kg
5	NMA10A89/T	SF	Titanium Four Hole Female Adapter, Axial (Offset by 7mm)	100kg
6	NMA10A60/A41	SF	Aluminium Double Pyramid Adapter 15mm	100kg
	NMA10A60/A50		Aluminium Double Pyramid Adapter 23mm	100kg
7	NTI-115DM	APC	Titanium Dual Pyramid Adapter 16mm	166kg
8	NMA11A16/T	SF	Titanium Pyramid Adapter with 30mm Core	125kg
	NMA11A16/S		Stainless Steel Pyramid Adapter with 30mm Core	150kg
9	NTA-1	TL	Titanium Rotation Adapter, 220g	150kg
			The NTA-1 rotation adapter can be assembled onto any 4-hole connector and can rotate up to 5°, due to its elongated mounting hole. An additional socket adapter is not required.	
10	NTA-2	TL	Titanium Rotation Adapter, 268g	150kg
			The NTA-2 rotation adapter offers even more adjustment options and up to 30° rotation with the M36 threaded connector.	
11	NMA12A7	SF	Socket Attachment Block - with poplar wood, 4 countersink head bolts 22A5/M6x20 including lamination dummy 12A3/E2 and sealing grease 22A100 with plastic for optimal sanding.	125kg
12	NAL-CPF-400	APC	Aluminium Round 4 Hole Laminating Plate Anodised BLACK (5/8" Thick with 1/2" Tie Groove) Includes Screws	136kg

Note: Whilst the components showing the blue water symbol are therefore waterproof, these components do not have drainage slots.

Manufacturer: SF= **Streifeneder** APC = **American Prosthetic Components** TL = **Trulife**