


## Structural and Socket Components - Waterproof with Drainage

**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.



### Titanium 30mm Structural Components

Part Number	Manufacturer	Description	Weight Limit	
1	NTI-10TC2-DH-SILVER	APC	262mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners.	166kg
2	NTI-17TC2-DH-SILVER	APC	440mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners.	166kg
3	NTI-10TC2-DH-BLACK	APC	262mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners, BLACK version.	166kg
4	NTI-17TC2-DH-BLACK	APC	440mm Aluminium Pylon with Titanium Receiver With drain holes and contains water resistant fasteners, BLACK version.	166kg
5	NTI-100C-DH	APC	Titanium Tube Clamp With drain holes and contains water resistant fasteners.	166kg
6	NTI-740L	APC	4 Prong Titanium Laminating Adapter with Pyramid	166kg
7	NTI-740-S	APC	4 Prong with Pryamid, Slots	166Kg

**Note:** The drainage slots included ensure these components are ideal for use on shower limbs, or for any application wherever there is submersion in, or frequent contact with water.

**Note:** See pages 252-253 for the Mix and Match APC Alignment Adapters Component Build Height Information

Manufacturer: APC = [American Prosthetic Components](#)