

Streifeneder Waterproof Shuttle-Lock

The Shuttle lock 15S2/W is characterised by its very small design and a high strength aluminium housing. Suitable for fresh, salt or even chlorinated water and manufactured in an excellent quality. The aluminium housing can be removed from the deep-drawing socket and used for the laminating socket. Processing is carried out using common deep-drawing or laminating methods.

Waterproof*



***Waterproof note:** Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Ordering Information

Waterproof Shuttle-Lock

| Part Number | Lock Type | Pin | Weight |
|-------------|--------------|----------|--------|
| 1 NMA15S2/W | Shuttle-Lock | Rastered | 90g |

Spare Parts Ordering Information

See page 280 for a full list of Series 150 spare parts

Features and Benefits

- No body weight restriction because the lock is not a weight bearing component
- Activity level: 1 to 4
- Can be used for thermoplastic and laminated sockets
- Aluminium housing
- Mounting height: approx 23mm
- Waterproof*
- Incl. Lamination dummy NMA10S1/E6
- Incl. special sealing grease NMA22A100
- Twist-proof pin guide made of stainless steel
- Locking unit completely exchangeable

