

Skinergy Plus provides an extremely high level of cosmesis for the below-knee user, incorporating a uniquely seamless finish and the latest innovation in silicone technology.


Because of its ease of use, increased durability, and improved realism, Skinergy Plus is a market leader for cosmetic off the shelf products.

Features and Benefits

- Integrated reinforcing mesh and advanced silicone material for a more durable finish
- Now available in 19 colour shades with unique TrueFinish™ micropigmentation and seamless finish for authentic realism
- Split or closed toe options to suit any foot size 220mm to 300mm
- Two calf sizes
- Easily removable for hardware adjustment, maintenance or repair
- Suitable for ankle disarticulation and partial foot applications
- Easy to clean using soap and water
- A Colour swatch kit is available to order to demonstrate all 19 shades
- 3-month warranty



Sandal and closed toe options



Semi-custom nail finish



Ordering Information

| Part Number | | Foot Size | Ankle Circ (cm) | | Mid Calf Circ (cm) | | Top of Cover Circ (cm) | | Approx Cover Length From Floor to: (cm) | |
|-------------|------------|-----------|-----------------|------|--------------------|------|------------------------|------|---|-----------|
| Split Toe | Closed Toe | Size | Min | Max | Min | Max | Min | Max | Mid Calf | Cover Top |
| TLC22LS/E* | TNS22LS/E* | 22-23 | 20 | 25 | 28 | 35 | 26 | 32.5 | 32 | 55 |
| TLC22RS/E* | TNS22RS/E* | 22-23 | 20 | 25 | 28 | 35 | 26 | 32.5 | 32 | 55 |
| TLC22LL/E* | TNS22LL/E* | 22-23 | 20 | 25 | 32 | 40 | 30 | 37.5 | 32 | 55 |
| TLC22RL/E* | TNS22RL/E* | 22-23 | 20 | 25 | 32 | 40 | 30 | 37.5 | 32 | 55 |
| TLC24LS/E* | TNS24LS/E* | 24-25 | 22 | 27.5 | 31 | 39 | 29 | 36 | 34.5 | 57 |
| TLC24RS/E* | TNS24RS/E* | 24-25 | 22 | 27.5 | 31 | 39 | 29 | 36 | 34.5 | 57 |
| TLC24LL/E* | TNS24LL/E* | 24-25 | 22 | 27.5 | 35 | 44 | 33 | 41 | 34.5 | 57 |
| TLC24RL/E* | TNS24RL/E* | 24-25 | 22 | 27.5 | 35 | 44 | 33 | 41 | 34.5 | 57 |
| TLC26LS/E* | TNS26LS/E* | 26-27 | 22 | 28 | 32 | 40 | 30 | 37.5 | 37 | 61 |
| TLC26RS/E* | TNS26RS/E* | 26-27 | 22 | 28 | 32 | 40 | 30 | 37.5 | 37 | 61 |
| TLC26LL/E* | TNS26LL/E* | 26-27 | 22 | 28 | 37 | 46 | 35 | 44 | 37 | 61 |
| TLC26RL/E* | TNS26RL/E* | 26-27 | 22 | 28 | 37 | 46 | 35 | 44 | 37 | 61 |
| TLC28LS/E* | TNS28LS/E* | 28-30 | 24 | 30 | 34 | 42.5 | 32 | 40 | 39.5 | 66 |
| TLC28RS/E* | TNS28RS/E* | 28-30 | 24 | 30 | 34 | 42.5 | 32 | 40 | 39.5 | 66 |
| TLC28LL/E* | TNS28LL/E* | 28-30 | 24 | 30 | 39 | 49 | 37 | 46 | 39.5 | 66 |
| TLC28RL/E* | TNS28RL/E* | 28-30 | 24 | 30 | 39 | 49 | 37 | 46 | 39.5 | 66 |

E* denotes colour. For example, TNS22LS/E5.

Warranty: 3 months

Colour Matching

| Part Number | Description |
|-------------|--------------------------------------|
| E24655 | Skinergy Plus True Finish Swatch Kit |

Skinergy covers can be requested with the following non-standard finishes: uncut, unpainted toenails. If a non-standard finish or bilateral cover is required, please specify when ordering.

TrueFinish Colours



Colour Matching:

These swatches are printed as guidance only. For colour matching, the swatch kits should be used (part number E24655).