

This lightweight polycentric knee is a robust unit with all the features you'd expect from a child's knee, but with one very useful additional feature. It has the facility to be set up as a semi auto or manual knee lock. This makes it particularly suitable for a child with bilateral limb loss, or with other functional issues, such as muscle weakness in the sound limb.

The lock can easily be removed should the child's ability not require it, or when an appropriate level of rehabilitation has been achieved. It can be just as easily replaced should there be a deterioration in function, such as following surgery.

Features and Benefits

- Polycentric design
- Lightweight
- Friction control of swing phase
- Extension assist
- Adjustable stance stability
- Optional manual lock
- Knee disarticulation adapter available (replaces the removable pyramid)



Ordering Information

Loaner units are available for this knee, please contact Customer Services for further details.

| Part Number | Description | Weight Limit |
|-------------|------------------------------------|--------------|
| NSSK610A | Child's Play Knee with Manual Lock | 55kg |

Additional Parts Ordering Information

| Image | Description | Weight Limit |
|-------|------------------------------------|--------------|
| | Lamination Bracket NSSK610A-DIS | 55kg |

| Image | Description | Weight Limit |
|-------|--|--------------|
| | Child's Play Energy Foot NSEF132, 144 | 65kg |

Additional Parts Ordering Information

| Image | Description | Weight Limit |
|-------|--|--------------|
| | Paediatric Rotatable Receiver - M36 Thread NJR-50FS-A | 60kg |

| Image | Description | Weight Limit |
|-------|--|--------------|
| | 3 Prong Laminating Base NSS-800BASE | 136kg |

Product Specification

| Weight Limit | Weight | Build Height | Warranty | Distance to Knee Centre |
|--------------|--------|--------------|----------|-------------------------|
| 55kg | 310g | 12cm | 2 years | 3cm |

Technical Information

Front View



Rear View

