


Paediatric Adapters - Structural and Socket Components



Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Socket Adapter with Rotation

Part Number	Manufacturer	Material	Weight	Weight Limit	
1	NMA100K4/A	SF	Aluminium	60g	45kg

Tube Clamp Adapter with Pyramid

Part Number	Manufacturer	Material	Weight	Weight Limit	
2	NMA110K15/A	SF	Aluminium	50g	45kg

Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Manufacturer: SF = **Streifeneder**

Double Adapter - Titanium

Part Number	Manufacturer	Material	Weight	Height	Weight Limit	
3	NJR-132D	APC	Titanium (pictured)	62g	32mm	60kg
	NJR-145D	APC	Titanium	74g	45mm	60kg
	NJR-160D	APC	Titanium	92g	60mm	60kg

Child's Rotatable Receiver

Part Number	Manufacturer	Material	Weight	Weight Limit	
4	NJR-50FS-A	APC	22mm Titanium Adapter	49g	60kg

Child's Lamination Plate

Part Number	Manufacturer	Material	Weight	Weight Limit	
5	NJR-AL-420L	APC	Aluminium Lamination Plate - supplied with M5 screws. Suitable for use with the Trulife adapter, NAAAK210.	51g	60kg

Child's 4-Hole Pyramid Adapter

Part Number	Manufacturer	Description	Weight	Weight Limit	
6	NJR-400B	APC	Titanium Paediatric 4-Hole Pyramid Adapter	36g	60kg
7	NJR-450BF-M6	APC	Titanium Paediatric 4-Hole Receiver Adapter	54g	60kg
8	NJR-400A	APC	Titanium Paediatric 4-Hole Pyramid, Axial Offset	45g	60kg
9	NJR-400D	APC	Titanium Paediatric 4-Hole Pyramid, Diagonal Offset	48g	60kg

Note: All paediatric 4 bolt components shown take M5 screws, but several of them are available in M6 screw versions. If required please speak to a Steeper Product Manager for further information.

Manufacturer: APC = **American Prosthetic Components**