


## Cosmetic Foam and Stockings - Structural and Socket Components



**Cosmetic Soft Foam Cover**  
For modular transtibial and knee disarticulation prostheses - made from PUR soft foam, applicable for left and right sides and with Ø 19mm bore.

Part Number	Manufacturer	Length
1 NMA600K7	SF	50cm

**Cosmetic Soft Foam Cover**  
For modular transfemoral and hip disarticulation prostheses - made from PUR soft foam with 30° knee flexion, a very distinctive pre-shaped calf, matching nylon cosmetic stocking (N200K15) and with Ø 19mm bore.

Part Number	Manufacturer	Side	Calf Circumference
2 NMA600K32/L	SF	Left	32cm
NMA600K32/R	SF	Right	32cm

**Cosmetic Stockings**  
For above-knee paediatric prostheses made from 100% polyamide (nylon).

Part Number	Manufacturer	Length	Foot Length	Prosthetic Foot Size
3 NMA200K15/1	SF	34cm	10cm	Up to 14
NMA200K15/2	SF	37cm	11cm	Up to 18
NMA200K15/3	SF	44cm	12cm	Up to 20

Manufacturer: SF = [Streifeneder](#)