


Paediatric Lamination Adapters - Structural and Socket Components



Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Child's Play Lamination Block - Titanium
If lamination is thick for strength, a longer screw may be required.

Part Number	Manufacturer	Description	Weight	Weight Limit
1 NSLB240	TL	Lamination Block	93g	65kg

Child's Play Foot Spacers

Part Number	Manufacturer	Foot Size	Height	Weight	Weight Limit
2 NSFS240	TL	13-15cm	6mm	6g	31kg
NSFS241	TL	13-15cm	25mm	20g	31kg
NSFS242	TL	16-18cm	6mm	9g	41kg
NSFS243	TL	16-18cm	25mm	28g	41kg
NSFS244	TL	19-21cm	6mm	11g	65kg
NSFS245	TL	19-21cm	30mm	45g	65kg

Child's Play Ankle Block

Part Number	Manufacturer	Foot Size	Weight	Weight Limit
3 NSAB314	TL	13-15cm	123g	35kg
NSAB315	TL	16-18cm	146g	45kg
NSAB316	TL	19-21cm	179g	65kg

Child's Play Symes Nut

Part Number	Manufacturer	Foot Size	Weight	Weight Limit
4 NSSY301	TL	13-15cm	14g	31kg
NSSY302	TL	16-18cm	20g	40kg
NSSY303	TL	19-21cm	28g	65kg

Child's Play Lamination Core

Part Number	Manufacturer	Foot Size	Weight	Weight Limit
5 NSLC301	TL	13-15cm	28g	31kg
NSLC302	TL	16-18cm	50g	40kg
NSLC303	TL	19-21cm	79g	65kg

Manufacturer: TL = Trulife