








## Paediatric Components - Structural and Socket Components



**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

### Paediatric Titanium College Park Structural Components

Part Number	Manufacturer	Description	Weight Limit	
1	N4HPR-TP	CP	4-Hole Pyramid Receiver, Standard	75kg
2	N4HPRR-TP	CP	4-Hole Pyramid Receiver, Rotatable	75kg
3	N4HPA-TP	CP	4-Hole Pyramid Adapter, Standard	75kg
4	N4HPAR-TP	CP	4-Hole Pyramid Adapter, Rotatable	75kg
5	NTA-T300-22	CP	Low-Profile Tube Adapter, 300mm length, 22mm diameter	60kg
6	NTC-T-22	CP	Tube Clamp, 22mm diameter	60kg
7	NTA-T300-25	CP	Low-Profile Tube Adapter, 300mm length, 25mm diameter	75kg
8	NTC-T-25	CP	Tube Clamp, 25mm diameter	75kg
9	N22K-PAP	CP	22mm Kit with Paediatric Pyramid Adapter (22mm Tube Adapter - 300mm length, 22mm Tube Clamp, 4-Hole Paediatric Pyramid Adapter)	60kg
10	N25K-PAP	CP	25mm Kit with Paediatric Pyramid Adapter (25mm Tube Adapter - 300mm length, 25mm Tube Clamp, 4-Hole Paediatric Pyramid Adapter)	75kg
11	N2530K-PAP	CP	25-30mm Kit with Paediatric Pyramid Adapter (25-30mm Tube Adapter - 300mm length, 25-30mm Tube Clamp, 4-Hole Paediatric Pyramid Adapter)	75kg
12	N4HLP-AP	CP	4-Hole Lamination Plate, Square, Aluminium	75kg

Manufacturer: CP = College Park

# Paediatric Components - Structural and Socket Components



**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Paediatric Aluminium 22mm Structural Components				
Part Number	Manufacturer	Description		Weight Limit
1	NJR-7 40L	APC	Paediatric 4 Prong Pyramid Adapter	60kg
2	NJR-10T2	APC	250mm Length 22mm Pylon	60kg
3	NJR-10TC2	APC	250mm Length 22mm Pylon with Titanium Receiver	60kg
3A	NJR-10TC2-I	APC	250mm Length 22mm Pylon with Short Adapter	60kg
4	NJR-100C	APC	22mm Titanium Adapter	60kg

**Note:** Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

Manufacturer: APC = American Prosthetic Components



Child's Tube Clamp with Adult Pyramid Receiver					
Part Number	Manufacturer	Material	Weight	Weight Limit	
5	NJR-100C-A	APC	22mm Titanium Adapter with 30mm Receiver	79g	60kg
6	NJR-10TC2-A	APC	Child's 22mm Tube with Adult Pyramid Receiver	-	60kg

Manufacturer: APC = American Prosthetic Components

## Paediatric Components - Structural and Socket Components



**Note:** Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

**Weight Limit Note:** Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Paediatric Titanium Pyramid and Receiver Adapters						
Part Number	Manufacturer	Description	Weight	Height	Weight Limit	
1	NJR-460-TRM	APC	Paediatric Clamping Base with Rotatable Pyramid	91g	13mm	60kg
2	NJR-465TRF	APC	Paediatric Clamping Base with Rotatable Receiver	94g	18mm	60kg
3	NJR-460TB	APC	Paediatric 4-Hole Threaded Base with Clamp	60g	9mm	60kg
4	NJR-00MS	APC	Paediatric Rotatable Pyramid	31g	5mm	60kg
5	NJR-50FS	APC	Paediatric Receiver, Female (pictured)	34g	10mm	60kg
	NJR-50FS-03	APC	Paediatric Receiver, Female	39g	13mm	60kg
	NJR-50FS-06	APC	Paediatric Receiver, Female	46g	16mm	60kg
	NJR-50FS-09	APC	Paediatric Receiver, Female	53g	19mm	60kg
	NJR-50FS-12	APC	Paediatric Receiver, Female	59g	22mm	60kg
6	NJR-500P-C	APC	Paediatric SACH Foot Adapter with Single, Centered Anti-Rotation Ridge. <b>Note:</b> Suitable for use with Trulife Child's Play Energy Foot	58g	7mm	60kg

**Note:** NJR-460-TRM, NJR-465TRF and NJR-460TB will bolt to NJR-AL-420L (page 336).

Manufacturer: APC = [American Prosthetic Components](#)