








Paediatric Components - Structural and Socket Components



Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Paediatric Titanium College Park Structural Components

| Part Number | Manufacturer | Description | Weight Limit | |
|-------------|--------------|-------------|--|------|
| 1 | N4HPR-TP | CP | 4-Hole Pyramid Receiver, Standard | 75kg |
| 2 | N4HPRR-TP | CP | 4-Hole Pyramid Receiver, Rotatable | 75kg |
| 3 | N4HPA-TP | CP | 4-Hole Pyramid Adapter, Standard | 75kg |
| 4 | N4HPAR-TP | CP | 4-Hole Pyramid Adapter, Rotatable | 75kg |
| 5 | NTA-T300-22 | CP | Low-Profile Tube Adapter, 300mm length, 22mm diameter | 60kg |
| 6 | NTC-T-22 | CP | Tube Clamp, 22mm diameter | 60kg |
| 7 | NTA-T300-25 | CP | Low-Profile Tube Adapter, 300mm length, 25mm diameter | 75kg |
| 8 | NTC-T-25 | CP | Tube Clamp, 25mm diameter | 75kg |
| 9 | N22K-PAP | CP | 22mm Kit with Paediatric Pyramid Adapter (22mm Tube Adapter - 300mm length, 22mm Tube Clamp, 4-Hole Paediatric Pyramid Adapter) | 60kg |
| 10 | N25K-PAP | CP | 25mm Kit with Paediatric Pyramid Adapter (25mm Tube Adapter - 300mm length, 25mm Tube Clamp, 4-Hole Paediatric Pyramid Adapter) | 75kg |
| 11 | N2530K-PAP | CP | 25-30mm Kit with Paediatric Pyramid Adapter (25-30mm Tube Adapter - 300mm length, 25-30mm Tube Clamp, 4-Hole Paediatric Pyramid Adapter) | 75kg |
| 12 | N4HLP-AP | CP | 4-Hole Lamination Plate, Square, Aluminium | 75kg |

Manufacturer: CP = College Park

Paediatric Components - Structural and Socket Components



Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

| Paediatric Aluminium 22mm Structural Components | | | | |
|---|--------------|-------------|--|--------------|
| Part Number | Manufacturer | Description | | Weight Limit |
| 1 | NJR-7 40L | APC | Paediatric 4 Prong Pyramid Adapter | 60kg |
| 2 | NJR-10T2 | APC | 250mm Length 22mm Pylon | 60kg |
| 3 | NJR-10TC2 | APC | 250mm Length 22mm Pylon with Titanium Receiver | 60kg |
| 3A | NJR-10TC2-I | APC | 250mm Length 22mm Pylon with Short Adapter | 60kg |
| 4 | NJR-100C | APC | 22mm Titanium Adapter | 60kg |

Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

Manufacturer: APC = American Prosthetic Components



| Child's Tube Clamp with Adult Pyramid Receiver | | | | | |
|--|--------------|----------|---|--------------|------|
| Part Number | Manufacturer | Material | Weight | Weight Limit | |
| 5 | NJR-100C-A | APC | 22mm Titanium Adapter with 30mm Receiver | 79g | 60kg |
| 6 | NJR-10TC2-A | APC | Child's 22mm Tube with Adult Pyramid Receiver | - | 60kg |

Manufacturer: APC = American Prosthetic Components

Paediatric Components - Structural and Socket Components



Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.


Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

| Paediatric Titanium Pyramid and Receiver Adapters | | | | | | |
|---|--------------|-------------|--|--------|--------------|------|
| Part Number | Manufacturer | Description | Weight | Height | Weight Limit | |
| 1 | NJR-460-TRM | APC | Paediatric Clamping Base with Rotatable Pyramid | 91g | 13mm | 60kg |
| 2 | NJR-465TRF | APC | Paediatric Clamping Base with Rotatable Receiver | 94g | 18mm | 60kg |
| 3 | NJR-460TB | APC | Paediatric 4-Hole Threaded Base with Clamp | 60g | 9mm | 60kg |
| 4 | NJR-00MS | APC | Paediatric Rotatable Pyramid | 31g | 5mm | 60kg |
| 5 | NJR-50FS | APC | Paediatric Receiver, Female (pictured) | 34g | 10mm | 60kg |
| | NJR-50FS-03 | APC | Paediatric Receiver, Female | 39g | 13mm | 60kg |
| | NJR-50FS-06 | APC | Paediatric Receiver, Female | 46g | 16mm | 60kg |
| | NJR-50FS-09 | APC | Paediatric Receiver, Female | 53g | 19mm | 60kg |
| | NJR-50FS-12 | APC | Paediatric Receiver, Female | 59g | 22mm | 60kg |
| 6 | NJR-500P-C | APC | Paediatric SACH Foot Adapter with Single, Centered Anti-Rotation Ridge. Note: Suitable for use with Trulife Child's Play Energy Foot | 58g | 7mm | 60kg |

Note: NJR-460-TRM, NJR-465TRF and NJR-460TB will bolt to NJR-AL-420L (page 336).

Manufacturer: APC = [American Prosthetic Components](#)

Paediatric Adapters - Structural and Socket Components



Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Socket Adapter with Rotation

| Part Number | Manufacturer | Material | Weight | Weight Limit | |
|-------------|--------------|----------|-----------|--------------|------|
| 1 | NMA100K4/A | SF | Aluminium | 60g | 45kg |

Tube Clamp Adapter with Pyramid

| Part Number | Manufacturer | Material | Weight | Weight Limit | |
|-------------|--------------|----------|-----------|--------------|------|
| 2 | NMA110K15/A | SF | Aluminium | 50g | 45kg |

Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Manufacturer: SF = **Streifeneder**

Double Adapter - Titanium

| Part Number | Manufacturer | Material | Weight | Height | Weight Limit | |
|-------------|--------------|----------|---------------------|--------|--------------|------|
| 3 | NJR-132D | APC | Titanium (pictured) | 62g | 32mm | 60kg |
| | NJR-145D | APC | Titanium | 74g | 45mm | 60kg |
| | NJR-160D | APC | Titanium | 92g | 60mm | 60kg |

Child's Rotatable Receiver

| Part Number | Manufacturer | Material | Weight | Weight Limit | |
|-------------|--------------|----------|-----------------------|--------------|------|
| 4 | NJR-50FS-A | APC | 22mm Titanium Adapter | 49g | 60kg |

Child's Lamination Plate

| Part Number | Manufacturer | Material | Weight | Weight Limit | |
|-------------|--------------|----------|--|--------------|------|
| 5 | NJR-AL-420L | APC | Aluminium Lamination Plate - supplied with M5 screws. Suitable for use with the Trulife adapter, NAAAK210. | 51g | 60kg |

Child's 4-Hole Pyramid Adapter

| Part Number | Manufacturer | Description | Weight | Weight Limit | |
|-------------|--------------|-------------|---|--------------|------|
| 6 | NJR-400B | APC | Titanium Paediatric 4-Hole Pyramid Adapter | 36g | 60kg |
| 7 | NJR-450BF-M6 | APC | Titanium Paediatric 4-Hole Receiver Adapter | 54g | 60kg |
| 8 | NJR-400A | APC | Titanium Paediatric 4-Hole Pyramid, Axial Offset | 45g | 60kg |
| 9 | NJR-400D | APC | Titanium Paediatric 4-Hole Pyramid, Diagonal Offset | 48g | 60kg |

Note: All paediatric 4 bolt components shown take M5 screws, but several of them are available in M6 screw versions. If required please speak to a Steeper Product Manager for further information.

Manufacturer: APC = **American Prosthetic Components**