

Breeze provides K2 users a cost-effective, yet highly-functional foot with a variety of features utilising College Park's high-end technology. The water-friendly system includes an optional water drainage port, corrosion-resistant materials and special modifications to clear water from the entire prosthetic system.

CES CP F 11: Use this code to read the Breeze™ clinical evaluation summary on the Steeper website.

Waterproof*

*Waterproof note: Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

- Low profile
- Intelliweave® composite for optimised movement
- Enviroshell foot shell design aids efficient toe roll-over
- Removable foot shell with anti-microbial properties, sandal toe as standard
- Water friendly with optional drainage port
- Integrated pyramid



Integrated pyramid

Optional water-drainage port



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

Low: Daily activities include mostly level ground walking, moving around the home and the community.

Ordering Example: NBE-LCEN-22-G3

Part ID	Side	Shell Colour	Mounting	Size	Firmness Category
NBE	L		EN		G 3
NBE	R		EN		G 3


Includes foot shell and CPI Spectra Sock.

Caucasian	C
Tan	T
Brown	B

21-30cm

Spare Parts Ordering Information

Foot Shell	Size (cm)
NBE-S ■ ● ▲	21-30



CAP	Size (cm)
NBE-C ■ ● -2125	21-25
NBE-C ■ ● -2630	26-30



- = "L" or "R" for Left or Right
- = Colour (C = Caucasian, T = Tan, B = Brown)
- ▲ = Foot size in cm (21-30)

Shelltread	Size (cm)
NBE-T ■ ● ▲	21-30



- = "L" or "R" for Left or Right
- = Colour (N = Neutral, J = Jet Black)
- ▲ = Foot size in cm (21-30)



CPI Spectra Sock	Size (cm)
NCPI-SS-SM	Small (19-25cm)
NCPI-SS-LG	Large (26-31cm)
NCPI-SSK-SM	Small 6 pack (19-25cm)
NCPI-SSK-LG	Large 6 pack (26-31cm)



Product Specification

Sizes	Mounting	Weight Limit	Build Height	Weight	Impact Level	Warranty	Foot Shell
21-25cm	Endo (IP)	100kg	6.2cm	523g Size 25cm including foot shell	L M H	2 years Foot shell - 3 months	C T B Caucasian, Tan or Brown Standard width
26-30cm		125kg	6.5cm				

Technical Information

