

Seattle SACH



The Seattle SACH foot is suitable for new amputees who are likely to ambulate with low to medium-low impact levels (K1). It comes fully assembled with an M10 pyramid and foot bolt as well as with a Lightfoot cosmesis in two colour options.

Features and Benefits

- Easy roll-over
- Lightweight
- Natural looking Lightfoot cosmesis
- Sandal toe
- Available in both light and dark foot cosmeses

Water-resistant*

*Water-resistant: Safe for contact with water and occasional submersion. Rinse immediately with fresh water and dry following exposure to water, particularly salt or chlorinated.

Cosmetics



Light/dark Sandal Toe Footshell

Ordering Information

Trial Period - There is no trial period on this foot.

Impact Level Descriptions

Low: Low speed walking and low intensity activities - no variable cadence.

Part Number	Description	Available Sizes
NSCH110	9.5mm Heel Light, Lightfoot cosmesis	22-29cm
NSCH113	9.5mm Heel Dark, Lightfoot cosmesis	22-29cm


Ordering Example - NSCH110-25L

Part Number	Size	Left/Right
NSCH110	25	L

Product Specification

Build Height and Heel Rise			
Part Number	Size (cm)	Build Height (cm)	Heel Rise (mm)
NSCH110	22-29	8.5	9.5
NSCH113	22-29	8.5	9.5

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
102kg	L M H	549g Size 25cm	22-29	9 Size 25-26 cm	2 years Foot shell - 1 year



Seattle SACH keel view



Seattle Natural



The Seattle Natural foot, with its Delrin Keel, has a slightly more responsive action than the SACH, making it ideal for low to medium-low impact levels, single speed, home ambulators or new amputees (K1-K2). It comes fully assembled with an M10 pyramid and foot bolt as well as a Lightfoot cosmesis, available in two colour options.

Features and Benefits

- Lightweight
- One standard Delrin Keel per foot size
- Propels the user through cadence
- Natural looking Lightfoot cosmesis
- Available in both Light and Dark foot cosmesis
- Sandal toe only

Water-resistant*

*Water-resistant: Safe for contact with water and occasional submersion. Rinse immediately with fresh water and dry following exposure to water, particularly salt or chlorinated.

Cosmetics



Light/dark Sandal Toe Footshell



Feet - Trulife Seattle SACH / Natural

Ordering Information

Trial Period - There is no trial period on this foot.

Impact Level Descriptions

Low: Low speed walking and low intensity activities - no variable cadence.

Part Number	Description	Available Sizes
NSNF160	9.5mm Heel Light, Lightfoot cosmesis	22-29 cm
NSNF163	9.5mm Heel Dark, Lightfoot cosmesis	22-29 cm


Ordering Example - NSNF160-25L

Part Number	Size	Left/Right
NSNF160	25	L

Product Specification

Build Height and Heel Rise			
Part Number	Size (cm)	Build Height (cm)	Heel Rise (mm)
NSNF160	22-29	8.5	9.5
NSNF163	22-29	8.5	9.5

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
102kg	L M H	448g Size 25cm	22-29	8.5	2 years Foot shell - 1 year



Seattle Natural keel view