

Sidekicks, are the first prosthetic feet designed specifically for use with foreshortened prostheses. Stride Control™ adjustment and replaceable bumpers allow dorsiflexion and plantarflexion in real-world environments, offering the ultimate in balance and stability for those with bilateral transfemoral amputation.

Sidekicks mimic anatomical ankle motion, allowing movement to occur naturally whilst reducing socket interface forces and eliminating unwanted torque. They're also gait-matched and approved for use in wet conditions.

CES CPI F12: Use this code to read the Sidekicks clinical evaluation summary on the Steeper website.

Waterproof*

*Waterproof note: Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

- Moulded tread is adhered to the foot with a molecular bond to maximise longevity
- Additional customisation possible through the use of Stride Control™ and replaceable bumpers
- Supported anatomical ankle motion and flexible platform for uneven ground compliance
- Waterproof*



Moulded tread



Mimics anatomical ankle rotation



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

Low: Daily activities include mostly level ground walking, moving around the home and the community.

Ordering Example: NSK-AK-2-G3

Part ID	Amp. Level	Quantity	Firmness Category
NSK	AK		G
		Single	1
		Pair	2


Firmness Category			
Weight (kg)	0-54	55-100	100-125
Size (cm)	One Size		
Low Impact	1	2	3

Spare Parts Ordering Information

Bumpers & Bushings	Description	Bumper Colour
NV-FBS-A	Front Bumper, Soft, Size A	S
NV-FBM-A	Front Bumper, Medium, Size A	M
NV-FBF-A	Front Bumper, Firm, Size A	F
NV-RBS-B	Rear Bumper, Soft, Size B	S
NV-RBM-B	Rear Bumper, Medium, Size B	M
NF-RBF-B	Rear Bumper, Firm, Size B	F
NV-RBX-B	Rear Bumper, Extra Firm, Size B	EF
NV-ABS-2530	Ankle Bushing, Soft, 25-30cm	S
NV-ABM-2530	Ankle Bushing, Medium, 25-30cm	M
NV-ABF-2530	Ankle Bushing, Firm, 25-30cm	F

Gait Match			
Firmness Category	Front Bumper	Ankle Bushings	Rear Bumper
1	S	S	S
2	M	M	M
3	F	F	F

S Soft M Medium F Firm EF Extra Firm



Product Specification

Mounting	Sizes	Weight Limit	Build Height
Endo (IP)	One Size	125kg	5.8cm
Weight of Unit		Impact	Warranty
422g		L M H	2 years

Technical Information

