

The Terrain foot was created to achieve anatomical ankle motion and energy return in a low-profile design. College Park use their tri-axial technology expertise in a new way. Versatile and service-free, the Terrain features a carbon fibre base and permanent bumpers inside the ankle housing for a progressive, comfortable response during low to high impact activities.

Waterproof*

*Waterproof note: Following exposure to chlorinated and salt water, rinse immediately with fresh water and dry.

Features and Benefits

- Tri-axial motion, a unique method achieving full movement on all three planes
- A more neutral feel due to its anatomical ankle position and pivot point
- Supported anatomical ankle motion and flexible platform for uneven ground compliance and a balanced energy return
- Maintenance free
- Waterproof*



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

| Ordering Example - NTN-LCEN-27-G6 | | | | | | Weight (kg) | Sizes (cm) | Low/Mod | High |
|-----------------------------------|------|--------------|----------|------|-------------------|-------------|------------|---------|------|
| Part ID | Side | Shell Colour | Mounting | Size | Firmness Category | | | | |
| NTN | L | | EN | | G | 0-65 | 21-30 | 1 | 2 |
| NTN | R | | EN | | G | 65-82 | | 2 | 3 |
| | | | | | | 82-100 | | 3 | 4 |
| | | | | | | 101-125 | | 4 | 5 |
| | | | | | | 126-150 | | 5 | 6 |
| | | | | | | 151-166 | 26-30 | 7 | 7 |

| | | |
|-----------|---|---------|
| Caucasian | C | 21-30cm |
| Tan | T | |
| Brown | B | |
| Jet Black | J | |

Spare Parts Ordering Information

| Foot Shell | Size (cm) |
|-------------|-----------|
| NHZ-S ■ ● ▲ | 21-30 |

| CAP | Size (cm) |
|-----------------|-----------|
| NBE-C ■ ● -2125 | 21-25 |
| NBE-C ■ ● -2630 | 26-30 |

■ = "L" or "R" for Left or Right
 ● = Colour (C = Caucasian, T = Tan, B = Brown, J = Jet Black)
 ▲ = Foot size in cm (21-30)

| Wedge Kit | Size (cm) |
|--------------|-----------|
| NHZ-HWK-2125 | 21-25 |
| NHZ-HWK-2630 | 26-30 |

| CPI Spectra Sock | Size (cm) |
|------------------|------------------------|
| NCPI-SS-SM | Small (19-25cm) |
| NCPI-SS-LG | Large (26-31cm) |
| NCPI-SSK-SM | Small 6 pack (19-25cm) |
| NCPI-SSK-LG | Large 6 pack (26-31cm) |

| Shelltread | Size (cm) |
|-------------|-----------|
| NHZ-T ■ N-▲ | 21-30 |
| Neutral | |

Product Specification

| Mounting | Sizes | Weight Limit | Build Height |
|-----------|---------|--------------|--------------|
| Endo (IP) | 21-30cm | 21-25cm | 149 kg |
| | | 26-30cm | 166 kg |
| | | | 8.6 cm |
| | | | 9 cm |

| Weight of Unit | Impact Level | Warranty | Foot Shell |
|----------------------------------|--------------|-------------------------------|---|
| 815 g Size 26 inc. foot shell | M H | 3 years foot shell - 6 months | C T B J Caucasian, Tan, Brown or Jet Black |
| | | | Standard Width |

