

With its excellent energy return and +/- 17 degrees of inversion and eversion, the Xtend Foot features characteristics comparable to a human foot. Constructed from tough glass fibre, carbon fibre and titanium this foot is ideal for all moderate-high activities including hiking, boating and golfing.

This foot provides a comfortable, natural walking style, even in challenging terrains like downward slopes and stairs. Suitable for users with above or below knee amputation and accommodates patients up to 150kg.



***Waterproof note:** Following exposure to chlorinated water, rinse immediately with fresh water and dry.

Features and Benefits

- Adapts to uneven surfaces through its lightweight, flexible structure
- +/- 17° of inversion and eversion
- Enables a more natural gait on slopes and inclines
- Strong and effective energy return in every step
- 150kg weight limit
- 3 year warranty (foot shell - 6 months)
- 30 day trial period



Waterproof*



Effective energy return



17° of Inversion and eversion

INTERACTIVE
Lindhe Xtend Foot Configurator



Scan To Configure



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions

Low: Daily activities include mostly level ground walking, moving around in the home and the community.

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

Weight Limit (kg)	Foot Size (cm)						
	23	24	25	26	27	28	29
60	NH100-323	NH100-324	-	-	-	-	-
80	NH100-423	NH100-424	NH100-425	NH100-426	NH100-427	NH100-428	NH100-429
100	NH100-523	NH100-524	NH100-525	NH100-526	NH100-527	NH100-528	NH100-529
125	-	-	NH100-625	NH100-626	NH100-627	NH100-628	NH100-629
150	-	-	-	NH100-726	NH100-727	NH100-728	NH100-729

How to Order:

For left or right, add **L** or **R** after the article code, and select footshell colour Caucasian, Brown, Jet Black or Tan (**FC, FB, FJ, or FT**) as a suffix. Example: **NH100-426L-FB**
All foot orders include Spectra Sock and footshell.

Replacement Foot Shells

For replacement footshells, replace **NH100** prefix with **NX900**. Example: **NX900-426L-FB**
To order please contact our Customer Services Team.

Product Specification

Material	Unit Weight	Weight Limit	Activity Level	Warranty
Glass Fibre, Carbon Fibre and Titanium	695g Size 26cm	150kg	Up to K3	Xtend Foot - 3 Years Foot Shell - 6 Months

Weight Limit (kg)	Foot Size (cm)																				
	23			24			25			26			27			28			29		
	A	B	C	A	B	C	A	B	C	A	B	C	A	B	C	A	B	C	A	B	C
60 kg	81	95	111	82	95	111	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
80 kg	81	95	111	82	96	112	80	104	120	90	106	122	88	106	122	89	109	125	92	109	125
100 kg	81	96	112	82	97	113	80	105	121	90	106	122	88	107	123	89	110	126	92	111	127
125 kg	-	-	-	-	-	-	80	106	122	90	107	123	88	107	123	89	112	128	92	112	128
150 kg	-	-	-	-	-	-	-	-	-	90	108	124	88	108	124	89	113	129	92	113	129

