

PROTEOR QUATTRO Package - Order Form

Included in the package:

- PROTEOR QUATTRO knee, including 5 or 6 year warranty* & the following accessories:
 - Wall Charger & Adapter
 - Battery Booster Kit
 - 1x foot: Exalto, Dynadapt, Highlander or Sierra
 - Set of adapters

*Warranty Period - NHS/Veteran = 5 Year Warranty Private = 6 Year Warranty

How to order a PROTEOR QUATTRO trial kit

Step 1.	Step 2.	Step 3.	Step 4 - Successful Trial Only
Submit a completed order form to customerservices@steepergroup.com Ensure a 'required by date' is highlighted if needed.	When the order form has been received, this will generate an order for a 30 day trial period for the QUATTRO MPK Package. You will not be charged or invoiced at this stage.	Following a successful trial, call Customer Services on 0870 240 4133 to obtain a RA number and return the full contents of the Trial QUATTRO knee box only . If the trial was unsuccessful, call Customer Services on 0870 240 4133 to obtain a RA number and return the full contents of the Trial QUATTRO knee box, Foot and Adapters .	After contacting Customer Services with your PO, the order will be finalised, delivered, and invoiced. Customer Services will confirm its estimated time of arrival.

N.B. If you do not require a trial and would like to place a definitive order, please submit this form, a purchase order using the QUATTRO part number and tick this box:

Patient Information:

Name: _____ Age: _____

Date knee is required for: _____ Weight (kg) : _____ Amputation Level (TF, KD, HD...): _____

Left: _____ Right: _____ Size (22-31cm): _____ Impact level**: Low: _____ Moderate: _____ High: _____ Bi-lateral _____

Notes: _____

Prosthetist: _____ Branch / Contractor: _____

**see opposite for impact level descriptions)

Adapters

The following adapters are included within the PROTEOR QUATTRO MPK package. Please select either Set A or Set B, for alternative adapters please place a separate order.

Set A (Formerly known as Set 1)

Description	Part Number	Weight Rating	Material
Tube Clamp, 30mm, Female Pyramid	NACC-00-12500-00	136kg	Titanium
Tube w/Titanium Female Pyramid, 30mm x 400mm	NACC-00-12800-00	136kg	Titanium
Pyramid, 4 Hole, Female	NACC-00-13200-00	136kg	Titanium

Set B (Formerly known as Set 2)

Description	Part Number	Weight Rating	Material
Tube Clamp, 30mm, Female Pyramid	NACC-00-12500-00	136kg	Titanium
Tube w/Titanium Female Pyramid, 30mm x 400mm	NACC-00-12800-00	136kg	Titanium
Base, 3 Prong, Rotating Pyramid	NACC-00-13800-00	136kg	Stainless Steel
Pyramid, Rotating, Female	NACC-00-14000-00	136kg	Titanium

Foot Selection - select one foot

Exalto (E) 



Available in sizes 22-31

Dynadapt (D) 



Available in sizes 23-30

Highlander (H) 



Available in sizes 23-30

Sierra (S) 



Available in sizes 23-30

Category Selection

Weight (kg)	Impact Level		
	Low	Mod	High
44-52	1	1	2
53-59	1	2	3
60-68	2	3	4
69-77	3	4	5
78-88	4	5	6
89-100	5	6	7
101-116	6	7	8
117-130	7	8	9
131-147	8	9	-
148-166	9	-	-



Impact Level Descriptions
Low: Daily activities include mostly ground-level walking, moving around in the home and the community.

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Size (cm)

22 23 24 25 26 27 28 29 30 31 Left: _____ Right: _____

Stiffness Category

1 2 3 4 5 6 7 8 9 Regular Toe: _____ Sandal Toe: _____

Skin Tone

Skin tone: Light: _____ Medium: _____ Dark: _____
*Sandal toe available in sizes 22-28 for Dynadapt, Highlander & Sierra Feet. Sandal toe is not available for the Exalto foot.

PROTEOR QUATTRO Warranty & Servicing

NHS/Veteran 5 years Private 6 years

Minimal servicing is required for QUATTRO, with no servicing needed until 36 months. For further warranty and servicing information please contact customer services.

QUATTRO Package

Please complete the form above, use the fields to specify your package part number and enter into the text box below:

Ordering Example - NQNX0010-NHS-EB

Part ID	Connection	Funding [†]	Foot Selection	Component Set
NQNX00				
Pyramid	10	NHS/Vet	NHS	Exalto E Set A A
Threaded	11	Private	6	Dynadapt D Set B B
				Highlander H
				Sierra S

[†]Note: NHS part numbers/ NHS Packages are only available to order through NHS or Veteran funding.

