

SoftskinAir® Range - Overview



**SoftSkinAir®  
NS30/NS40/NS50**  
Breathable Silicone Liners



**SoftSkinAir®  
NS30CL/NS40CL/NS50CL**  
Breathable Silicone Liners



**SoftSkinAir®  
2L**  
Breathable Silicone Liners



**SoftSkinAir®  
SPORT**  
Breathable Silicone Liner



Padding  
Stability

**SoftSkinAir®  
Cushion 30**



Padding  
Stability

**SoftSkinAir®  
Cushion 40**



**SoftSkinAir®  
TF**  
Comfortable & Breathable  
Transfemoral Liners

Material	Silicone (with micropores)				Silicone (with micropores)		Silicone (with micropores)		Silicone (with micropores)		Silicone (with micropores)	
Shore Hardness	30/40/50		30/40/50		30/50		40		30		40	
Thickness	3mm		3mm		3mm		3mm		3mm		3mm	
Length	370mm		370mm		370mm		370mm		370mm (23.5 - 28mm) 400mm (28 - 36mm) 430mm (36mm<)		350mm	
Amputation Level	Transtibial		Transtibial/Transfemoral		Transtibial		Transtibial		Transtibial/Transfemoral		Transfemoral	
Distal Connection	Yes		No		Yes		Yes		No		Yes	
Mobility Grades	NS30 & NS40 1,2 & 3	NS50 2,3 & 4	NS30CL & NS40CL 1,2 & 3	NS50CL 2,3 & 4	1,2 & 3		2,3 & 4		1,2 & 3		1,2,3 & 4	
Warranty	6 months		6 months		6 months		6 months		6 months		6 months	