


The Streifeneder ContexGel.sleeves ensure correct suspension and fit in lower and upper limb prosthetic applications. They're effective at keeping air out of elevated system sockets and keeping the prosthesis firmly in place. Smaller sizes are suitable for upper limb below-elbow prosthetic sockets; larger sizes are suitable for below-knee amputees.

The latest bi-elastic knit design prevents tension and pinching on the knee joint, with 3mm antibacterial Polymer Gel for optimal residual limb health.

Features and Benefits

- Conical shaped sleeve
- Available in two colours: mixed anthracite or sand
- Round-knit design with no lateral seams
- Soft, abrasion-resistant textile cover which can be easily donned or doffed
- Features an outer matrix to eliminate elongation in the distal and proximal area
- Excellent circumferential elasticity compensates for residual limb volume changes



Antibacterial Technology (SteriOne)



Ordering Information

Colour: Mixed Anthracite				
Size	Distal	Proximal	Length	Part Number
S	15 - 23cm	23 - 35cm	37cm	NMA30L10/S
M	21 - 31cm	29 - 44cm	37cm	NMA30L10/M
L	28 - 43cm	36 - 52cm	38.5cm	NMA30L10/L
XL	34 - 48cm	42 - 58cm	38.5cm	NMA30L10/XL
XXL	38 - 52cm	43 - 64cm	38.5cm	NMA30L10/XXL
XXXL	42 - 56cm	50 - 70cm	38.5cm	NMA30L10/XXXL

Colour: Sand				
Size	Distal	Proximal	Length	Part Number
S	15 - 23cm	23 - 35cm	37cm	NMA30L50/S
M	21 - 31cm	29 - 44cm	37cm	NMA30L50/M
L	28 - 43cm	36 - 52cm	38.5cm	NMA30L50/L
XL	34 - 48cm	42 - 58cm	38.5cm	NMA30L50/XL
XXL	38 - 52cm	43 - 64cm	38.5cm	NMA30L50/XXL
XXXL	42 - 56cm	50 - 70cm	38.5cm	NMA30L10/XXXL

ContexGel.protect



To avoid premature tearing of the ContexGel knee sleeves, we recommend using the ContexGel.protect as a protective cover for sharp-edged prosthetic sockets. This protective cover is made of durable and expandable fabric and has an inner silicone edge on the distal side to prevent slipping.

ContexGel.protect cover				
Size	Side	For Knee Sleeve	Length	Part Number
S	Universal	Size S	19cm	NMA3S50/S
M	Universal	Size M	19cm	NMA3S50/M
L	Universal	Size L	22cm	NMA3S50/L
XL	Universal	Size XL & XXL	22cm	NMA3S50/XL

Size Information

ContexGel Knee Sleeves are available in six sizes. When determining patient size, please position the sleeve so that the middle section is sat over the knee/elbow joint.

To determine the optimal size of the knee sleeve, measure the distal socket circumference and the proximal thigh circumference (approx. 15cm above the centre of the knee joint).

Protective Inner Sleeve



These sleeve protectors from Össur are made of durable fabric that extends the life of the suspension sleeve and protects it from damage caused by the socket's edge. Silicone pimples on the bottom edge of the sleeve hold it securely in place on the socket.

Protective Sleeve	
Size	Part Number
Small	NM000202
Medium	NM000203
Large	NM000204

Technical Information

Material	Polymer Gel
Length	37 - 38.5cm
Amputation Level	Transtibial
Warranty	3 months