


SureFit Sleeve

The SureFit Sleeve is constructed using High-Density TPE Gel (HD Gel) to provide maximum comfort for the patient and long-lasting security. With a 30° pre-flexion it allows for ease when bending and reduces bunching behind the knee; whilst the intelligent weave is adapted to differing areas of stretch and wear.

Features and Benefits

- 30° pre-flexion
- 3mm of uniform thickness
- Low modulus of elasticity for correct compression and stability
- One length (40m) to accommodate most users in three sizes
- Two colour options available: Grey or Beige



Innovative weave design



High-Density Gel





Sleeves - Steeper - SureFit Sleeve

Ordering Information

Circumference	Sleeve Size*	Colour	Product Part No.
20 - 34cm	Small	Beige	NSFX20-SM-B
20 - 34cm	Small	Grey	NSFX20-SM-G
31 - 40cm	Medium	Beige	NSFX26-MD-B
31 - 40cm	Medium	Grey	NSFX26-MD-G
35 - 60cm	Large	Beige	NSFX32-LG-B
35 - 60cm	Large	Grey	NSFX32-LG-G

Size Information

*Measure the thigh circumference 10cm above the patella. If the measured circumference falls within two sleeve sizes, please order the smaller size.



Technical Information

Material	High-Density Thermoplastic Elastomer (TPE) Gel
Length	40cm
Warranty	6 months

