

This polyester and Lycra® residual limb sock utilises COOLMAX® fibres to wick away perspiration for all-day comfort.

Features and Benefits

- Available in 3 and 5 ply
- Soft and lightweight material ensures user comfort and reduces irritation
- Non-chafing and breathable characteristics, ideal for active users
- Available with a distal hole.
- Machine wash and air dry (do not use laundry detergent)

Ordering Information

1. Please order in multiple of six
2. Determine part number using the chart below
3. To order with a distal hole, replace the 'N5' prefix with 'NH' - e.g. **NH1135U12**

Application Recommendation

Ideal for users who have delicate/sensitive residual limb skin, or require perspiration control.

Ordering Example

3 ply Cool Sock II, 7" Top Width, 12" Length = **N51135U12**



Length	Top Width (Proximal)							
	15-18 cm		18-20 cm		20-23 cm		22-25 cm	
	Toe Width (Distal)							
	7-10 cm	10-13 cm	10-13 cm	13-15 cm	15-18 cm			
	5 ply	5 ply	3 ply	5 ply	3 ply	5 ply	5 ply	
20-25 cm	N51153U08	N51154U08	N51135U08	N51155U08	N51136U08	N51156U08	N51157U08	
25-31 cm	N51153U10	N51154U10	N51135U10	N51155U10	N51136U10	N51156U10	N51157U10	
31-36 cm	N51153U12	N51154U12	N51135U12	N51155U12	N51136U12	N51156U12	N51157U12	
40-46cm	N51153U16	N51154U16	N51135U16	N51155U16	N51136U16	N51156U16	N51157U16	
51-56 cm	N51153U20	N51154U20	N51135U20	N51155U20	N51136U20	N51156U20	N51157U20	
61-66 cm	N51153U24	N51154U24	N51135U24	N51155U24	N51136U24	N51156U24	N51157U24	

Ply Identification Colour Code

Yellow	3 ply
Green	5 ply