

# Transfemoral Compression Socks

Intelligent and effective products with deep functional effects. The Össur compression seal identifies selected high-tech products from Össur in medicine, lifestyle and sports.


The manufacturing process provides accurately dosed compression on the body - for positive effects on health, performance and well-being.

## Features and Benefits

- Soft, smooth material with a comfortable feel for all-day wearing
- Movable AK waistbelt



Clima Comfort




Clima Fresh



# Socks - Össur® - Transfemoral Compression Sock

## Size Selection

### Above-knee



## Ordering Information

Length	Size	Cir. cF (cm)	Cir. cG (cm)	Part No.
Short 18-22cm	1	30-32	42-45	N1882 ▲ 11
	2	32-35	45-49	N1882 ▲ 12
	3	35-38	49-53	N1882 ▲ 13
	4	38-41	53-57	N1882 ▲ 14
	5	41-44	57-61	N1882 ▲ 15
	6	44-47	61-65	N1882 ▲ 16
Medium 23-27cm	1	30-32	42-45	N1882 ▲ 21
	2	32-35	45-49	N1882 ▲ 22
	3	35-38	49-53	N1882 ▲ 23
	4	38-41	53-57	N1882 ▲ 24
	5	41-44	57-61	N1882 ▲ 25
	6	44-47	61-65	N1882 ▲ 26
Long 28-32cm	1	30-32	42-45	N1882 ▲ 31
	2	32-35	45-49	N1882 ▲ 32
	3	35-38	49-53	N1882 ▲ 33
	4	38-41	53-57	N1882 ▲ 34
	5	41-44	57-61	N1882 ▲ 35
	6	44-47	61-65	N1882 ▲ 36
X-Long 33-37cm	1	30-32	42-45	N1882 ▲ 41
	2	32-35	45-49	N1882 ▲ 42
	3	35-38	49-53	N1882 ▲ 43
	4	38-41	53-57	N1882 ▲ 44
	5	41-44	57-61	N1882 ▲ 45
	6	44-47	61-65	N1882 ▲ 46

▲ - Denotes Compression, add 1 for CCL1 and 2 for CCL2

## Technical Information

### Indications

- Reduction of oedema/swelling of the lower limb after amputation
- Maintenance of limb volume when the prosthesis is removed
- For fixing soft tissues

### Contraindications

- If the medical findings forbid compression (e.g. heart failure)

Material	Umbrellan® with Clima-Comfort®
Size / Length	See Ordering Information
Amputation Level	Knee disarticulation / transfemoral
Cleaning	30° hand wash
Compression Class	CCL1 and CCL2