

Transtibial Compression Socks

Intelligent and effective products with deep functional effects for optimum comfort. The Össur compression seal identifies selected high-tech products from Össur in medicine, lifestyle and sports.

The manufacturing process ensures that the body receives precisely dosed compression, which has a positive impact on health, performance, and well-being.

Features and Benefits


- Soft, smooth material for comfortable all-day wearing
- Silicone beaded BK top-band for superior suspension
- Smaller BK size



Socks - Össur® - Transtibial Compression Socks

Size Selection

Below-knee



Technical Information

Indications

- Reduction of oedema/swelling of the lower limb after amputation
- Maintenance of limb volume when the prosthesis is removed
- For fixing soft tissues

Contraindications

- If the medical findings forbid compression (e.g. heart failure)

Material	Umbrellan® with Clima-Comfort®
Size / Length	See Ordering Information
Amputation Level	Transtibial
Cleaning	30° hand wash
Compression Class	CCL1 and CCL2

Ordering Information

Length	Size	Cir. cC (cm)	Cir. cE (cm)	Cir. cF (cm)	Part No.
Short 28-34cm	0	26-28	28-30	38-40	N1881 ▲10
	1	28-30	30-32	40-42	N1881 ▲11
	2	30-33	32-35	42-45	N1881 ▲12
	3	33-36	35-38	45-48	N1881 ▲13
	4	36-39	38-41	48-51	N1881 ▲14
	5	39-42	41-44	51-54	N1881 ▲15
Medium 35-41cm	6	42-45	44-47	54-57	N1881 ▲16
	0	26-28	28-30	38-40	N1881 ▲20
	1	28-30	30-32	40-42	N1881 ▲21
	2	30-33	32-35	42-45	N1881 ▲22
	3	33-36	35-38	45-48	N1881 ▲23
	4	36-39	38-41	48-51	N1881 ▲24
Long 42-48cm	5	39-42	41-44	51-54	N1881 ▲25
	6	42-45	44-47	54-57	N1881 ▲26
	0	26-28	28-30	38-40	N1881 ▲30
	1	28-30	30-32	40-42	N1881 ▲31
	2	30-33	32-35	42-45	N1881 ▲32
	3	33-36	35-38	45-48	N1881 ▲33
	4	36-39	38-41	48-51	N1881 ▲34
	5	39-42	41-44	51-54	N1881 ▲35
	6	42-45	44-47	54-57	N1881 ▲36

▲ - Denotes compression, add 1 for CCL1 and 2 for CCL2