

## Espire Classic Plus

The Espire Classic Plus is a mechanical elbow featuring a manual lock and forearm counterbalance mechanism. Using a plug-in board at the top of the arm, prosthetists can attach various inputs and batteries negating the need for external cabling. These plug-ins are routed through the elbow to the distal end for attachment of a wider range of terminal devices. This model provides a discreet compartment for housing the S-Charge System or other manufacturers' batteries without having to modify the arm.

### Technical Specifications

**WEIGHT**

Small | 665g  
Standard | 685g

**DIAL PLACEMENT OPTIONS**

Medial  
Lateral

**SIDES**

Left  
Right

**WARRANTY**

1 year



## Espire Classic

The Espire Classic is a mechanical elbow that features a forearm counterbalance mechanism and manual lock. Contains no electronics or circuitry to provide a lightweight option.

### Technical Specifications

**WEIGHT**

Small | 651g  
Standard | 671g

**DIAL PLACEMENT OPTIONS**

Medial  
Lateral

**SIDES**

Left  
Right

**WARRANTY**

1 year

## Espire Basic

The Espire Basic is a body-powered elbow, featuring a manual lock and an economical dial-adjustable spring-assist to enhance ease of movement. Contains no electronics or circuitry to provide a lightweight device.

### Technical Specifications

**WEIGHT**

Small | 649g  
Standard | 669g

**DIAL PLACEMENT OPTIONS**

Medial  
Lateral

**SIDES**

Left  
Right

**WARRANTY**

1 year



# Espire™ Elbow

Classic Plus, Classic, and Basic

With an emphasis on function and design, the Espire Elbow range is easy to use and anatomically correct. Three of the five models are mechanical joints with varying features. The Espire Classic Plus, Classic, and Basic are sleek, lightweight, and intuitive solutions.

## Features

### Electronic Pass-through | Classic Plus

To allow for myoelectric control of terminal devices, the Classic Plus has a plug-in board at the top of the arm where prosthetists can attach various batteries and inputs without the need for external cabling. Compatible with the Steeper S-Charge System.



### Anatomical Shape | Classic Plus, Classic, and Basic

This design blends technology with humanity to appeal to the wide range of patient demographics. The Espire is designed to emulate an anatomical forearm and wrist in every way possible. Because of this, the arm looks natural in a comfortable resting position.



### Comfortable Counterbalance | Classic Plus and Classic

The Classic and Classic Plus offer this adjustable feature that counterbalances the weight of the forearm, wrist, and terminal device, negating the effects of gravity. This allows the patient to position the forearm with minimal effort.



	Number of inputs/ outputs	Locking functions	Elbow function	Custom shell options	Clean-line shell and anatomical shape	Wrist options
<b>Classic Plus</b>	4/4	Manual	Counterbalance	●	●	Small (45mm) & Standard (50mm)
<b>Classic</b>	0		Counterbalance	●	●	
<b>Basic</b>	0		Spring-assist	●	●	

### Colour Options

Custom colours also available



Unit 3 Stourton Link, Intermezzo Drive,  
Leeds, LS10 1DF  
Email: [customerservices@steepergroup.com](mailto:customerservices@steepergroup.com)  
Tel: +44 (0) 870 240 4133

[www.steepergroup.com](http://www.steepergroup.com)  
STPPRI64 Issue 1 December 2020  
© 2020 Steeper Group. All rights reserved.

 **steeper**  
Creating life's turning points, together